



China Dental Lab Invisible Strong U Shape Lingual Wire Retainer Single Strand Comfort Level for Orthodontic Treatment

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians




Product Specification

- Usage: Orthodontic Treatment
- Suitable For: All Ages
- Durability: Strong
- Removability: Non-removable
- Bonding Method: Direct Bonding
- Wire Configuration: Single Strand
- Aesthetics: Invisible
- Comfort Level: High
- Highlight: **U shape lingual wire retainer, single strand orthodontic retainer, invisible dental lab retainer**



More Images



Product Description

Product Description:

Our Pr

One of the unique features of this retainer is its high level of comfort. Unlike traditional retainers that can be bulky and uncomfortable, the Bonded Lingual Wire Retainer is designed to be as discreet and comfortable as possible. The wire is carefully contoured to fit the lingual surface of the teeth, ensuring that patients can speak and eat normally while wearing the retainer.

Another key feature of the Bonded Lingual Wire Retainer is its aesthetics. Because the wire is bonded to the lingual surface of the teeth, it is virtually invisible. This makes it an excellent option for patients who want to maintain a natural-looking smile while ensuring their teeth remain in their newly aligned positions.

Overall, the Bonded Lingual Wire Retainer is an excellent choice for patients who want a strong, durable, and comfortable retainer that is virtually invisible. Whether you are a teenager or an adult, this lingual bonded retainer is an excellent option for maintaining a healthy, beautiful smile for years to come.

Features:

Product Name: Lingual Wire Retainer(U shape)

Wire Type: Round

Usage: Orthodontic Treatment

Aesthetics: Invisible

Material: Stainless Steel

Suitable For: All Ages

Keywords: Lingual Retainer Bonded, Bonded Lingual Wire, Bonded Lingual Wire

Technical Parameters:

Lingual Wire Retainer	
Wire Configuration	Single Strand
Usage	Orthodontic Treatment
Durability	Strong
Removability	Non-removable
Wire Type	Round
Material	Stainless Steel
Aesthetics	Invisible
Bonding Method	Direct Bonding
Suitable For	All Ages
Comfort Level	High

Applications:

The ADL Lingual Wire Retainer is a high-quality orthodontic product that is suitable for all ages. Made of durable stainless steel, this lingual bonded retainer is designed for direct bonding and is available in a round wire type. It has been certified by ISO, CE, and FDA, ensuring its safety and efficacy.

The Lingual Retention Wire is commonly used in orthodontic treatment where it is attached to the lingual side of the teeth. It is designed to maintain the position of teeth after treatment, preventing them from shifting back to their original position.

The ADL Lingual Wire Retainer is ideal for patients who have undergone orthodontic treatment and want to maintain their newly aligned teeth. It is also suitable for patients who have had braces and are at risk of their teeth shifting back to their original position.

The Lingual Retention Wire is commonly used in scenarios where patients want to maintain their newly aligned teeth without the use of visible retainers. It is also commonly used in scenarios where patients have a high risk of their teeth shifting back, such as those with a history of teeth grinding or jaw clenching.

The ADL Lingual Wire Retainer has a minimum order quantity that is negotiable, and the price is also negotiable. It has a delivery time of 72 hours and accepts payment terms of Paypal or T/T. With a supply ability of 630 technicians, the Lingual Retention Wire is readily available for all patients.

In conclusion, the ADL Lingual Wire Retainer is an effective and convenient solution for maintaining newly aligned teeth, and it is suitable for patients of all ages. With its durable stainless steel material and direct bonding method, it is a reliable and safe option for orthodontic treatment.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China