

## Non Removable Lingual Wire Retainer Invisible Direct Bonded Lingual Retainer Wire

Our Product Introduction

for more products please visit us on [adidentalab.com](http://adidentalab.com)

### Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Material: Stainless Steel
- Suitable For: All Ages
- Bonding Method: Direct Bonding
- Durability: Strong
- Wire Type: Round
- Comfort Level: High
- Usage: Orthodontic Treatment
- Wire Configuration: Single Strand
- Highlight: Non Removable Lingual Wire Retainer, Lingual Wire Retainer Invisible, Direct Bonded Lingual Retainer Wire



### More Images



### Product Description

Non-removable Lingual Wire Retainer Invisible Direct Bonding

**Product Description:**

Our Product

The Bonded Lingual Wire is made of high-quality materials that are designed to withstand the wear and tear of daily use. Whether you are eating, speaking, or simply going about your daily routine, this retainer will remain firmly in place, providing constant support for your teeth. This makes it an ideal choice for patients who are at risk of relapse or who have experienced teeth shifting in the past. In addition to its strong and durable design, the Lingual Wire Retainer is also a highly effective orthodontic treatment option. By providing constant pressure on the teeth, the retainer helps to maintain their alignment over time, reducing the risk of relapse and ensuring that your smile stays straight and healthy for years to come. Overall, the Lingual Wire Retainer is an excellent choice for anyone looking for a strong and durable orthodontic treatment option that is suitable for patients of all ages. Whether you are a teenager or an adult, this retainer can help you achieve the beautiful, healthy smile you deserve.

### Features:

**Product Name:** Lingual Wire Retainer

**Material:** Stainless Steel

**Suitable For:** All Ages

**Aesthetics:** Invisible

**Bonding Method:** Direct Bonding

**Wire Type:** Round

This Lingual Wire Retainer is a great option for those looking for an invisible and convenient solution. Made from durable stainless steel, it is suitable for all ages and uses a direct bonding method for easy application. The round wire type ensures a comfortable fit, making the Lingual Retainer Bonded a great option for long-term use.

### Technical Parameters:

Durability	Strong
Usage	Orthodontic Treatment
Wire Configuration	Single Strand
Removability	Non-removable
Material	Stainless Steel
Suitable For	All Ages
Bonding Method	Direct Bonding
Comfort Level	High
Aesthetics	Invisible
Wire Type	Round

### Applications:

The Bonded Lingual Wire Retainer is a highly effective and comfortable solution for maintaining your teeth's alignment. This product is bonded to the back of your teeth, making it virtually invisible and providing an aesthetically pleasing option for orthodontic treatment. The Lingual Retainer Bonded is suitable for people of all ages who have undergone orthodontic treatment and are looking for a comfortable and effective way to maintain their teeth's alignment.

The ADL Lingual Wire Retainer is an excellent choice for those who want to maintain their teeth's alignment and improve their overall oral health. This product is certified by ISO/CE/FDA, ensuring that it meets the highest quality standards. The minimum order quantity for this product is negotiable, and the price is negotiable as well. The delivery time for this product is 72 hours, and payment terms include PayPal or T/T.

The ADL Lingual Wire Retainer is made possible by our team of 630 technicians who are dedicated to providing high-quality orthodontic products. The Lingual Wire Retainer is incredibly comfortable, making it an excellent choice for those who want to maintain their teeth's alignment without sacrificing comfort. The bonding method used for this product is direct bonding, ensuring a secure and long-lasting bond.

The Lingual Wire Retainer is perfect for a wide range of occasions and scenarios. Whether you are looking to maintain your teeth's alignment after orthodontic treatment, or you are looking to improve your overall oral health, the ADL Lingual Wire Retainer is an excellent choice. This product is suitable for all ages and is incredibly comfortable, making it an ideal choice for anyone who wants to maintain their beautiful smile.



**Advanced Dental Laboratory**



+86 18665847934



[info@adldentalab.com](mailto:info@adldentalab.com)



[adldentalab.com](http://adldentalab.com)

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China