



Titanium Locator Dental Attachments For Overdentures In Silver And Gold Finish

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Cad/Cam Compatibility: Yes
- Material: Titanium
- Color: Silver And Gold
- Easy To Clean: Yes
- Usage: Useful Attachment For Overdentures
- Durability: Long-lasting
- Highlight: **titanium locator dental attachments, titanium overdenture locator attachments, overdentures locator dental attachments**



Product Description

Titanium Locator Attachment for Overdentures in Silver and Gold Finish

Product Description:

One of the best things about the Locator Attachment product is that it is incredibly easy to clean. This makes it a great product for anyone who wants to maintain good oral hygiene. The silver and gold color of the Locator Extension Attachment also adds a touch of elegance to your overdenture.

Another great thing about the Locator Attachment product is that it provides a secure fit. This is thanks to the Locator Extension Attachment, which keeps your overdenture in place. You can trust that your overdenture will not shift or move around when you use the Locator Attachment product.

The Locator Attachment product is also very easy to use. The Locator Connection Kit makes it easy to attach the Locator Extension Attachment to your overdenture. You don't need any special tools or equipment to use this product.

Overall, the Locator Attachment product is an excellent investment for anyone who uses overdentures. It is durable, easy to clean, and provides a secure fit. The Locator Extension Attachment is made of titanium, which means that it will last for a long time.

Features:

Product Name: Locator Attachment

Cad/Cam Compatibility: Yes

Durability: Long-lasting

Color: Silver And Gold

Usage: Useful Attachment For Overdentures

Easy to Clean: Yes

Technical Parameters:

Product Attribute	Description
Usage	Useful Attachment For Overdentures
Easy to Clean	Yes
Color	Silver And Gold
Durability	Long-lasting
Material	Titanium
Cad/Cam Compatibility	Yes
Product Name	Locator Extension Attachment
Alternate Names	Attachment for Locators, Attachment for Locators

Applications:

The Locator Mounting Attachment comes in plastic bags and box packaging, and has a delivery time of 72 hours. The payment terms for this product include Paypal or T/T, and we have a supply ability of 630 technicians.

This attachment is made of titanium and comes in silver and gold colors. It is compatible with Cad/Cam and is very easy to clean. The product is a useful attachment for overdentures as it provides a great connection between the denture and the implant.

The Locator Mounting Attachment is a great product for dental clinics, hospitals, and laboratories. It can be used in various occasions such as dental implant surgeries, overdenture fittings, and dental prosthetics.

The Locator Connection Kit is a great accessory for the Locator Mounting Attachment. It contains all the necessary components for a successful attachment of the Locator onto the overdenture.

The Attachment for Locators is a great way to secure overdentures in place. It is highly recommended for patients who struggle with loose dentures. The Locator Mounting Attachment provides great stability and comfort to the patient.



Advanced Dental Laboratory



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China