



China Dental Lab Adjustable Patient Specific Keyway Precision Attachment with PFM

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Ease Of Use: Easy To Insert And Remove
- Customization: Can Be Customized For Individual Patient Needs
- Color: Silver
- Retentive Force: Adjustable
- Size: Various Sizes Available
- Hygiene: Easy To Clean And Maintain
- Compatibility: Fits Most Dental Implants
- Durability: Long-lasting
- Highlight: dental lab precision attachment, adjustable patient specific keyway, PFM precision dental attachment



More Images



Product Description

Product Description from China dental lab Advanced Dental Lab:

Our F

One of the standout features of the Key&Key precision attachment is its retention mechanism. The key and keyway design ensures a secure fit for the partial denture, preventing any unwanted movement or discomfort for the patient.

Another advantage of the Key&Key Way Precision Attachment Partial Denture is its hygiene. The attachment is easy to clean and maintain, making it a convenient and practical choice for both patients and dental professionals alike. But perhaps the most impressive aspect of the Key&Key Way Precision Attachment is its customization. Each attachment can be tailored to meet the individual needs of the patient, ensuring a comfortable and effective fit. This level of personalization is essential for patients who require partial dentures, as it can greatly improve their quality of life. Overall, the Key&Key Way Precision Attachment Partial Denture is an excellent choice for anyone in need of a reliable and customizable precision attachment. Whether you're a dental professional looking to provide the best possible care for your patients, or a patient in need of a comfortable and effective partial denture, the Key&Key precision attachment is an ideal solution.

Features:

Product Name: Key&Key Way Precision Attachment

Retentive Force: Adjustable

Type: Precision Attachment

Size: Various Sizes Available

Color: Silver

Ease Of Use: Easy To Insert And Remove

This Key&Key Way Precision Attachment is a precision attachment device that provides adjustable retentive force for partial dentures with keyway attachment. It comes in various sizes and a silver color, and is easy to insert and remove.

Technical Parameters:

Product Attribute	Description
Product Name	Key&Key Precision Attachment
Type	Precision Attachment
Size	Various Sizes Available
Hygiene	Easy To Clean And Maintain
Compatibility	Fits Most Dental Implants
Precision	High Precision
Color	Silver
Retentive Force	Adjustable
Customization	Can Be Customized For Individual Patient Needs
Ease Of Use	Easy To Insert And Remove
Versatility	Can Be Used For Multiple Dental Prostheses

Applications:

One of the primary applications of the Key&Key Way Precision Attachment Partial Denture is for patients who have lost one or more teeth and want to restore their smile. This product is perfect for those who want a long-lasting, comfortable solution that is easy to insert and remove.

Another scenario where the Key&Key Way Precision Attachment Partial Denture is useful is for patients who require partial dentures but do not want to use traditional metal clasps. The Key&Key Way Precision Attachment Partial Denture is a great alternative, as it is more aesthetically pleasing and easier to maintain.

The Key&Key Way Precision Attachment Partial Denture is also ideal for patients who have sensitive gums or suffer from periodontal disease. These patients require a high level of hygiene, and the Key&Key Way Precision Attachment Partial Denture is easy to clean and maintain, ensuring good oral health.

The ADL Key&Key Way Precision Attachment Partial Denture comes in plastic bags and boxes, with a delivery time of 72 hours. The minimum order quantity, price, and payment terms are negotiable, making it accessible to a wide range of dental practices. The product is supplied by 630 technicians who are trained to ensure that it meets the highest standards of quality.

In conclusion, the Key&Key Way Precision Attachment Partial Denture is an excellent product for anyone who wants a comfortable, durable, and aesthetically pleasing solution for missing teeth. Its adjustable retentive force, long-lasting durability, and easy-to-use design make it a top choice for dental practices around the world.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China