



## Telescope Inner Co-Cr And Outer Co-Cr Metal Partial With Composite

Our Product Introduction

for more products please visit us on [adldentalab.com](http://adldentalab.com)

### Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Biocompatibility: Yes
- Material: Inner Co-Cr And Outer Co-Cr Metal Partial With Composite
- Color: White
- Chemical Resistance: Excellent
- Durability: High
- Application: Dental/Orthodontic
- Highlight: **Co-Cr Metal Partial, Telescope Metal Partial, Composite Metal Partial**



### More Images



## Product Description

### Telescope Inner Co-Cr and Outer Co-Cr Metal Partial with Composite

#### Product Description:

The Telescope is designed specifically for dental and orthodontic applications. It is an excellent choice for patients who require a partial denture but want a prosthesis that is durable and long-lasting. It is easy to take out and clean. This prosthesis is made from a combination of Inner Co-Cr and Outer Co-Cr Metal Partial with Composite, which are both materials that are known for their long lasting.

One of the standout features of the Telescope is its high level of durability. This prosthesis is designed to withstand the wear and tear of everyday use, making it an excellent choice for patients who need a prosthesis that will last for many years. It is also resistant to staining and discoloration, which helps to keep it looking bright and white for longer.

Overall, the Telescope is an excellent choice for patients who require a telescope prosthesis that is durable, long-lasting, biocompatible, and aesthetically pleasing. Whether you are a dental professional or a patient, the Telescope is a high-quality option that will meet all of your needs.

#### Features:

**Product Name:** Telescope

**Corrosion Resistance:** High

**Durability:** High

**Application:** Dental

**Chemical Resistance:** Excellent

**Can be used as a Telescope partial denture**

**Can be used as a Telescope prosthesis**

**Also suitable for use as a Telescope partial denture**

#### Technical Parameters:

| Technical Parameter | Description  |
|---------------------|--|
| Material            | Inner Co-Cr and Outer Co-Cr Metal Partial with Composite |
| Size                | Customizable   |
| Durability          | High   |
| Biocompatibility    | Yes  |
| Application         | Dental   |
| Weight              | Lightweight  |
| Chemical Resistance | Excellent  |
| Color               | White  |

#### Applications:

ADL Telescope (Inner Co-Cr and Outer Co-Cr Metal Partial with Composite) is a perfect fit for dental and orthodontic applications as it is highly resistant to chemicals. The product is available in plastic bags and boxes, and its delivery time is only 72 hours. Payment for the product can be made through Paypal or T/T, and the minimum order quantity and price are negotiable.

The telescope dental prosthesis is an application that requires a high level of precision and durability. ADL Telescope has been designed to meet these requirements, making it the perfect solution for dental and orthodontic applications.

The product's lightweight nature also makes it easy to use and handle, ensuring that the dentist can perform their work with ease and comfort.

ADL Telescope is supplied by 630 technicians, ensuring that the product is always available for use. The product's packaging ensures that it is protected during transportation and arrives at the destination safely.



**Advanced Dental Laboratory**

+86 18665847934

info@adldentalab.com

adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China

