



Convenient Occlusal Sport Guard Bite Splint Design Easy To Clean Comfortable Protective Device For Teeth Grinding

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Durability: Durable
- Easy To Clean: Yes
- Comfortable: Yes
- Adjustability: Adjustable For Fit
- Convenient: Yes
- Flexibility: High
- Design: Bite Splint Design
- Breathability: Breathable
- Highlight: occlusal splint for teeth grinding,
sport guard bite splint,
comfortable teeth grinding protector



More Images



Product Description

Product Description:

Our F

The Sport Guard is an advanced dental appliance designed to provide superior protection and comfort for individuals who require effective oral care during physical activities. Crafted with precision using cutting-edge technology, this product stands out for its exceptional quality, durability, and user-friendly features. Sport guard is used to shield your teeth during high-impact activities, this product offers an ideal solution that combines convenience, effectiveness, and comfort.

One of the most notable features of this Sport Guard is its material composition. This durability makes the guard a cost-effective investment, as it maintains its integrity and performance over an extended period without frequent replacements.

Adjustability is a critical aspect when it comes to oral appliances, and the Sport Guard excels in this regard. The product is designed to be adjustable for fit, allowing users to customize the guard to their unique dental structure. This adjustability ensures a secure and comfortable fit, reducing any irritation or discomfort that might arise from using a poorly fitting guard.

Effectiveness is paramount in any dental guard, and this Sport Guard delivers exceptional performance. For athletes, it offers robust protection against impacts that could otherwise cause dental injuries.

Convenience is also a key benefit of this product. The ease of adjusting the fit, combined with its durable and breathable design, makes the Sport Guard an appliance that fits seamlessly into daily routines. Users can quickly adapt the guard to their needs without requiring professional assistance.

In summary, Sport Guard is an innovative and reliable dental appliance that leverages the advantages of sport guard.

Features:

Breathable design ensures comfort during wear

High flexibility for ease of movement and comfort

Adjustable for fit to provide personalized comfort

Offer precise and custom fit

Lightweight

Technical Parameters:

Easy To Clean	Yes
Convenient	Yes
Flexibility	High
Effective	Yes
Breathability	Breathable
Comfortable	Yes
Durability	Durable
Design	Bite Splint Design
Adjustability	Adjustable For Fit
Usage	Protection during sport

Applications:

The ADL Sport Guard is expertly crafted in Shenzhen, China, combining advanced technology with high-quality materials to deliver superior protection and comfort. Featuring a design, this product offers exceptional durability and flexibility, making it an ideal choice for various application occasions and scenarios.

In sports and athletic activities, the ADL Sport Guard plays a crucial role in protecting athletes from dental injuries. Its high flexibility and durable construction absorb impact effectively, safeguarding teeth and gums during contact sports such as football, basketball, and martial arts. The breathable material allows for comfortable wear even during intense physical exertion, while the bite splint design ensures that the guard stays securely in place.

The product's design is certified by ISO, CE, and FDA, ensuring it meets international standards for safety and quality.

With a flexible supply ability of 630 technicians and a negotiable minimum order quantity and price, ADL offers customized solutions tailored to individual needs. The product is conveniently packaged in plastic bags and boxes, with a swift delivery time of 72 hours. Payment can be made securely via Paypal or T/T, making the purchasing process smooth and efficient.

In summary, the ADL Sport Guard is a versatile and reliable solution for dental protection and sports safety. Its breathable, durable, and easy-to-clean design, ensures optimal performance in a wide range of occasions and scenarios, from clinical dental applications to active sports environments.

FAQ:

Q1: What is the brand name of the Sport Guard?

A1: The Sport Guard is branded as ADL.

Q2: Where is the Sport Guard manufactured?

A2: It is manufactured in Shenzhen, China.

Q3: What certifications does the ADL Sport Guard have?

A3: The product is certified with ISO, CE, and FDA approvals.

Q4: What is the minimum order quantity for the ADL Sport Guard?

A4: The minimum order quantity is negotiable, so you can discuss your needs with the supplier.

Q5: How long does it take to deliver the Sport Guard after placing an order?

A5: The delivery time is typically 72 hours after the order is confirmed.

Q6: What are the payment terms for purchasing the ADL Sport Guard?

A6: Payment can be made via Paypal or T/T (Telegraphic Transfer).

Q7: How is the Sport Guard packaged?

A7: The product is packaged in plastic bags and boxes to ensure safe delivery.

Q8: What is the supply ability of the ADL Sport Guard?

A8: The supplier has a supply ability of 630 technicians to meet production demands efficiently.



Advanced Dental Laboratory



+86 18665847934



info@adlentalab.com



adlentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China