



China Dental Lab Titanium And Zirconia All on X in 16 shades including A1 A2 A3 from dental lab for natural dental aesthetics

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians




Product Specification

- Turnaround Time: In 72 Hours
- Design: All On X Design
- Materials Used: Titanium And Zirconia
- Stain Resistance: High
- Shade: A1/A2/A3.....16 Shade
- Product Type: Dental Laboratory Services
- Customization: CAD/CAM Designed For Precision Fit
- Shipping Options: UPS/DHL/FEDEX



More Images



Product Description

Our F

Product Description:

When it comes to high-quality dental restorations, partnering with a reliable outsource dental lab is essential for dental professionals seeking precision, durability, and aesthetic excellence. Our China dental lab-Advanced Dental Lab stands out as a premier choice, offering cutting-edge dental solutions that cater to a wide range of prosthetic needs. Specializing in advanced designs such as the All On X system, including the popular All On 4 technique, our lab ensures that dental practitioners receive products that meet the highest standards of quality and functionality.

One of the primary technical parameters that set our China dental lab-Advanced Dental Lab apart is the use of premium materials such as high standard titanium and zirconia. Titanium is renowned for its exceptional strength, biocompatibility, and corrosion resistance, making it an ideal choice for implant frameworks and overdentures. Zirconia, on the other hand, offers outstanding aesthetics, durability, and a natural tooth-like appearance, perfect for fixed prostheses. By combining these materials, our lab produces restorations that are not only long-lasting but also visually appealing, ensuring patient satisfaction and confidence.

Our lab is designed to support a broad range of dental applications, including both overdentures and fixed prostheses. This versatility allows dental professionals to offer personalized treatment plans tailored to each patient's unique dental needs. Whether you require a sturdy overdenture solution or a fixed prosthetic designed for optimal function and appearance, our China dental lab-Advanced Dental Lab has the expertise and technology to deliver superior results. The All On X design approach, which encompasses various implant placement strategies including All On 4, serves as a foundation for many of our restorations, providing stability and excellent load distribution for patients.

Efficiency and rapid turnaround times are critical in the dental industry, where timely delivery can significantly impact patient care and satisfaction. Understanding this necessity, our China dental lab guarantees a turnaround time of just 72 hours. This swift processing time ensures that dental professionals can maintain their treatment schedules without unnecessary delays, enhancing patient experience and practice efficiency. Our streamlined workflow and advanced manufacturing technologies enable us to meet tight deadlines while maintaining uncompromising quality standards.

Choosing to outsource dental lab work to our China dental lab-Advanced Dental Lab means leveraging our extensive technical expertise and state-of-the-art manufacturing capabilities. Our team of skilled technicians employs the latest digital design tools and production methods to create precise, custom-fit restorations. The All On X design methodology, including the widely adopted All On 4 implant concept, is expertly implemented to optimize implant placement and prosthetic stability, resulting in improved functional outcomes and patient comfort.

In summary, our China dental lab-Advanced Dental Lab offers a comprehensive solution for dental professionals looking to outsource dental lab services without sacrificing quality or turnaround time. With a focus on advanced materials like titanium and zirconia, support for both overdentures and fixed prostheses, and a commitment to rapid 72-hour delivery, we provide a reliable and efficient partner for dental restorations. The incorporation of the All On X design, including the All On 4 system, underscores our dedication to modern, effective dental implant solutions that meet the evolving needs of patients and clinicians alike.

For dental practices aiming to enhance their restorative offerings, our China dental lab-Advanced Dental Lab represents a strategic choice. By outsourcing dental lab work to us, you gain access to cutting-edge designs, premium materials, and expedited service, all crafted to ensure optimal patient outcomes and satisfaction. Trust our expertise to deliver excellence in every prosthesis, helping you elevate your practice and provide outstanding care.

Applications:

The ADL brand, a renowned China dental lab with over 13 years of experience in the dental lab industry, offers exceptional products designed for various dental application occasions and scenarios. As a trusted outsource dental lab, ADL specializes in providing high-quality dental solutions, including the widely popular all on 6 dental implants. With certification from ISO, CE, and FDA, ADL ensures that every product meets stringent international standards for safety and performance.

One of the primary application occasions for ADL's products is in comprehensive dental restoration cases, particularly for patients requiring full-arch rehabilitation. The all on 6 dental implant solution is ideal for scenarios where patients need a strong, stable, and long-lasting prosthesis. Utilizing milled titanium bar connections, the products boast a bending strength of $\geq 1000\text{Mpa}$, providing superior durability and reliability in clinical applications. This makes ADL's offerings perfect for both permanent and temporary restorations in demanding dental environments.

Dental clinics and practitioners looking to outsource dental lab work can greatly benefit from partnering with ADL. The China dental lab's supply ability of 630 skilled technicians ensures timely and efficient production, with a guaranteed turnaround time of just 72 hours. This quick delivery is crucial for cases that require fast processing without compromising quality. Additionally, flexible minimum order quantities and negotiable pricing accommodate both small and large-scale dental practices, making ADL a versatile partner for diverse dental projects.

Packaging is handled with care, using plastic bags and boxes to protect the integrity of each dental component during shipping. The payment terms are convenient, accepting Paypal or T/T, facilitating smooth transactions for international clients. ADL's high-strength dental products, backed by extensive industry experience and rigorous quality control, cater to a wide range of dental restoration scenarios, from single implants to complex all on 6 solutions.

In summary, ADL, the leading China dental lab, is the optimal choice for dental professionals seeking a reliable outsource dental lab capable of delivering high-quality, certified dental products with fast turnaround times. Whether it's for routine dental restorations or advanced all on 6 implant cases, ADL's expertise, robust product strength, and efficient service make it an indispensable partner in modern dentistry.

Customization:

ADL is a leading outsource dental lab based in China, specializing in high-quality dental products with ISO, CE, and FDA certifications. Our China Dental Lab offers customizable sizes to support overdentures and fixed prostheses, including all on 4 solutions, ensuring a high-strength and precise fit. With a supply ability of 630 skilled technicians, we guarantee a fast turnaround time with delivery within 72 hours.

We provide a variety of shades from A1 to A16 to meet diverse aesthetic needs, and our products are CAD/CAM designed for precision fit. Packaging is done carefully using plastic bags and boxes to ensure product safety. Our minimum order quantity is negotiable, and pricing is flexible to accommodate your specific requirements.

Payment terms include Paypal or T/T for your convenience. Choose ADL China Dental Lab for reliable, high-quality dental products with excellent customization services tailored to your needs.

Packing and Shipping:

Our China dental lab-Advanced Dental Lab products are carefully packaged to ensure maximum protection during transit. Each item is securely wrapped with cushioning materials and placed in a sturdy, custom-designed box to prevent damage. The packaging is labeled clearly with product details and handling instructions to facilitate smooth delivery.

For shipping, we offer reliable and timely options including air freight, sea freight, and express courier services. We collaborate with trusted logistics partners to guarantee that your dental lab products arrive safely and on schedule. Tracking information is provided once the order is dispatched, allowing you to monitor the shipment until it reaches your doorstep.



Advanced Dental Laboratory



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China