



China Dental Lab Compatible With Most Dental Implant Systems All on 4 Offering In 72 Hours Turnaround Time And Array Technical Parameters

Our Product Introduction

for more products please visit us on addentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Size: Customizable
- Experience: Over 13 Years In Dental Lab Industry
- Technical Parameters: Array
- Stain Resistance: High
- Connection Type: Milled Titanium Bar Connection
- Product Type: Dental Laboratory Services
- Design: All On X Design
- Bending Strength: $\geq 1000\text{Mpa}$
- Highlight: Dental Implant Systems, China Dental Lab, 72 Hours Turnaround Time Dental Lab



More Images



Product Description

Our F

Product Description:

China Dental Lab stands out as one of the top 5 dental labs in the industry, offering unparalleled expertise and quality in dental prosthetics and restorations. With over 13 years of experience in the dental lab industry, this lab has established a strong reputation for delivering outstanding products that combine durability, precision, and aesthetic appeal. Their comprehensive understanding of dental materials and technologies ensures that every product meets the highest standards demanded by dental professionals worldwide. One of the key attributes of China Dental Lab's products is their impressive bending strength, which is ≥ 1000 MPa. This remarkable strength guarantees that dental restorations such as crowns, bridges, and frameworks are not only resilient but also built to withstand the rigors of daily use. Patients benefit from restorations that are both long-lasting and reliable, contributing to improved oral health and satisfaction. The high bending strength also makes these products suitable for complex procedures, including the popular all on x implant solutions, where durability and stability are critical.

In addition to mechanical strength, China Dental Lab emphasizes the importance of aesthetic quality through high stain resistance. Their dental restorations maintain their natural appearance over time, resisting discoloration from food, beverages, and other common staining agents. This feature ensures that patients enjoy a bright, natural-looking smile for years, which is essential for confidence and comfort. The combination of high stain resistance with expert craftsmanship results in products that not only perform well but also look exceptional.

Customization is another area where China Dental Lab excels. Utilizing state-of-the-art CAD/CAM technology, every restoration is designed for a precise fit tailored to each patient's unique dental anatomy. This advanced digital design process minimizes errors and adjustments, leading to faster turnaround time and enhanced patient satisfaction. Dentists can rely on this precision fit to streamline their workflow and provide superior care. Whether producing single crowns or full arch restorations, the lab's CAD/CAM designed products ensure optimal functionality and comfort.

Understanding the critical nature of timely delivery in the dental field, China Dental Lab offers reliable shipping options through major carriers such as UPS, DHL, and FEDEX. This ensures that dental professionals receive their orders promptly, enabling them to meet patient needs without delay. The lab's commitment to fast turnaround time means that even complex cases can be completed efficiently, supporting smooth clinical operations and patient treatment schedules.

China Dental Lab's dedication to quality, innovation, and customer service has firmly positioned it among the top 5 dental labs globally. Their comprehensive product offerings, combined with advanced manufacturing techniques and robust quality control, make them a trusted partner for dental professionals seeking excellence. From high-strength materials to precise CAD/CAM customization and dependable shipping, every aspect of their service is designed to support superior dental care.

In summary, China Dental Lab offers a compelling combination of high bending strength (≥ 1000 MPa), extensive industry experience exceeding 13 years, and superior stain resistance. Their CAD/CAM designed restorations guarantee precision fit, while flexible shipping options through UPS, DHL, and FEDEX ensure fast turnaround time. Trusted by dental professionals worldwide for all on x and other advanced restorative procedures, China Dental Lab remains a leader in delivering durable, aesthetic, and precisely crafted dental solutions.

Applications:

The ADL dental lab, a top 5 dental lab based in China, offers exceptional products and services tailored for a variety of dental application occasions and scenarios. Renowned for its high-quality dental solutions, ADL is certified by ISO, CE, and FDA, ensuring compliance with international standards. This makes it an ideal choice for dental professionals seeking a reliable outsource dental lab to support their practice needs.

One of the primary application occasions for ADL products is the All On 6 dental implant procedure. ADL specializes in the All On X design, providing durable and precise dental prosthetics with a bending strength of ≥ 1000 Mpa. This strength is crucial for long-lasting restorations that can withstand the stresses of daily use. The availability of 16 shades, including A1, A2, and A3, allows for highly aesthetic and customized dental restorations that match patients' natural teeth perfectly.

ADL's supply ability of 630 technicians ensures timely production and delivery, with a standard delivery time of just 72 hours. This quick turnaround is vital for dental clinics looking to outsource dental lab work efficiently without compromising quality. The packaging details, including plastic bags and boxes, ensure safe transportation, while shipping options via UPS, DHL, and FEDEX provide flexibility and reliability in global logistics.

Dental practitioners can negotiate minimum order quantities and pricing, which is competitively structured with volume discounts to support both small and large-scale operations. Payment terms are convenient, accepting PayPal or T/T, making transactions smooth and secure. The combination of competitive pricing, rapid delivery, and certified quality makes ADL a preferred partner for dental labs and clinics worldwide.

In summary, ADL's China-based dental lab is perfectly suited for various dental application scenarios, especially for clinics performing All On 6 implant restorations. As a top 5 dental lab option, ADL provides comprehensive outsource dental lab services that combine advanced design, aesthetic customization, and robust mechanical properties, meeting the evolving demands of modern dentistry.

Customization:

ADL is a top 5 dental lab based in China, specializing in high-quality dental products with certifications including ISO, CE, and FDA. Our China dental lab offers comprehensive product customization services tailored to meet your specific needs.

We use premium materials such as Titanium and Zirconia to ensure durability and biocompatibility. Our designs focus on the All On X technique, providing excellent compatibility with most dental implant systems.

With a supply ability supported by 630 skilled technicians, we guarantee high stain resistance and superior product quality. Packaging is done carefully using plastic bags and boxes to ensure safe delivery.

Our minimum order quantity is negotiable, with pricing also flexible to accommodate your requirements. We offer fast delivery within 72 hours and convenient payment terms through Paypal or T/T.

Shipping options include reliable carriers such as UPS, DHL, and FEDEX, ensuring your customized dental products arrive promptly and securely.

Packing and Shipping:

Our China Dental Lab products are carefully packaged to ensure maximum protection during transit. Each item is securely placed in customized, sturdy packaging materials that prevent damage and contamination. The packaging is designed to maintain the integrity and hygiene of the dental products, complying with international safety and quality standards.

For shipping, we offer reliable and efficient logistics services with options for air, sea, and express delivery. All shipments are tracked and insured to guarantee safe and timely arrival at your destination. Our experienced shipping team coordinates closely with customers to provide seamless delivery and excellent after-sales support.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong

