



## China Dental Lab High Esthetics Customized Full Acrylic Denture Full Service for Worldwide

Our Product Introduction

for more products please visit us on [addentalab.com](http://addentalab.com)

### Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Turnaround Time: In 72 Hours
- Feature: Removable Denture
- Materials Used: BPS Acrylic/Yamahachi Teeth
- Product Type: China Dental Laboratory Services
- Shade: 16 Shades
- Shipping Options: UPS/DHL/FEDEX
- Experience: Over 13 Years In Dental Lab Industry
- Pricingmodel: Competitive Pricing With Volume Discounts
- Advantages: Precise And Cost Effective
- Highlight: China Dental Lab Full Acrylic Denture, Customized Full Acrylic Denture, High Esthetics Full Acrylic Denture



### More Images



### Product Description

**Product Description:**

Our F

China Dental Lab ---Advanced Dental Lab, top 5 dental lab in China, is a leading provider in the dental industry, offering high-quality dental solutions with over 13 years of extensive experience in the dental lab sector. Our commitment to excellence and innovation has established us as a trusted partner for dental professionals worldwide. Specializing in removable dentures, our products are designed to deliver superior comfort, durability, and esthetics, ensuring patient satisfaction at every stage.

One of the standout features of our dental lab services is compatibility with high esthetics requirements. We understand that aesthetics play a crucial role in dental restorations, especially when it comes to removable dentures. Our products are meticulously crafted to match natural teeth, providing a seamless and natural appearance. With a wide range of shades available, including A1, A2, A3, and extending to a comprehensive palette of 16 shades, we ensure that each denture perfectly complements the patient's unique dental profile. Our expertise is not limited to conventional dentures; we are also proficient in providing solutions for advanced implant techniques such as all on 6 and all on x systems. These cutting-edge implant protocols require precision and customization, and our dental lab is equipped to meet these demanding specifications. Whether it's all on 6 or all on x implants, our removable dentures are designed to integrate flawlessly with implant frameworks, offering both functionality and esthetic appeal.

At China Dental Lab, we pride ourselves on offering competitive pricing without compromising quality. Our pricing model is structured to provide affordability while supporting dental practices of all sizes. We offer volume discounts to ensure that clinics and dental professionals can access premium dental lab services at reduced costs when ordering in bulk. This approach not only fosters long-term partnerships but also empowers dental providers to deliver the best care to their patients without budget constraints.

Our extensive experience of over 13 years in the dental lab industry has enabled us to refine our processes and adopt the latest technologies. This experience translates into superior craftsmanship, reliable turnaround times, and consistent quality control. We stay updated with global dental trends and standards, allowing us to support advanced treatment plans such as all on 6 and all on x implant restorations with confidence and precision.

The removable denture products from China Dental Lab are designed with patient comfort and functionality in mind. Using high-quality materials and advanced manufacturing techniques, we ensure that our dentures are lightweight, durable, and easy to maintain. Our removable dentures provide excellent fit and retention, enhancing the patient's ability to speak, chew, and smile naturally. Furthermore, our extensive shade range from A1 to a full 16-shade spectrum allows dental professionals to customize the esthetic outcome to match the patient's natural dentition perfectly.

In summary, China Dental Lab stands out as a premier dental laboratory offering specialized removable dentures compatible with high esthetic demands and advanced implant protocols such as all on 6 and all on x. With over 13 years of industry experience, competitive pricing structures including volume discounts, and a comprehensive shade range, we are dedicated to supporting dental professionals worldwide. By choosing China Dental Lab, you partner with a reliable, innovative, and cost-effective dental lab committed to excellence in every restoration we produce.

### Applications:

The ADL dental lab product, proudly made in China, is designed to meet the diverse needs of dental professionals and patients alike. With over 13 years of experience in the dental lab industry, ADL has established itself as a trusted brand known for quality and reliability. Certified by ISO, CE, and FDA, ADL ensures that every product meets stringent international standards, making it a preferred choice for dental labs worldwide.

One of the standout features of ADL's dental products is their compatibility with high esthetics, particularly suitable for creating removable dentures using premium materials such as BPS Acrylic and VITA teeth. This combination guarantees a natural look and durable performance, ensuring patient satisfaction. Available in a wide range of shades from A1 to A16, these dentures can be customized to perfectly match the natural teeth of patients, enhancing their smile and confidence.

ADL's dental lab products are particularly advantageous for scenarios requiring a fast turnaround time. With a delivery time of just 72 hours, dental professionals can rely on ADL to provide timely solutions, which is crucial in urgent cases or busy clinics. The minimum order quantity is negotiable, providing flexibility for both small and large scale dental labs. Packaging is thoughtfully done using plastic bags and boxes to ensure the safety and hygiene of the products during transit.

The product is highly applicable in various dental treatment occasions, including the popular "all on 6" implant technique, where precision and esthetics are paramount. ADL's skilled team of 630 technicians ensures meticulous craftsmanship for each case, supporting dental practitioners in delivering excellent results. Payment terms are convenient with options like Paypal or T/T, facilitating smooth transactions globally.

In summary, ADL from China offers a comprehensive dental lab solution that caters to a broad range of application occasions and scenarios, from routine removable dentures to advanced implant-supported restorations like all on 6. Its fast turnaround time, quality certifications, customizable shade options, and experienced technical team make ADL an ideal partner for dental professionals seeking reliable and esthetically superior dental products.

### Support and Services:

Our China Dental Lab product is backed by comprehensive technical support and services to ensure optimal performance and customer satisfaction. We provide expert guidance on installation, operation, and maintenance to help you maximize the product's capabilities. Our skilled support team is available to assist with troubleshooting and resolving any technical issues promptly. Additionally, we offer regular software updates and training sessions to keep your lab up-to-date with the latest advancements. Trust our dedicated service to enhance your dental lab's efficiency and productivity.

### Packing and Shipping:

Our China Dental Lab products are carefully packaged to ensure maximum protection during transit. Each item is securely wrapped with high-quality materials and placed in sturdy, custom-designed boxes to prevent any damage.

We use shock-absorbing and moisture-resistant packaging to maintain the integrity and hygiene of dental products throughout the shipping process.

All shipments are tracked with reliable courier services, offering timely delivery and real-time updates to our customers worldwide.

For bulk orders, we provide reinforced pallet packaging to guarantee safe handling and efficient transportation.

Our commitment is to deliver your dental lab products in perfect condition, ready for immediate use upon arrival.



**Advanced Dental Laboratory**



+86 18665847934



[info@adldentalab.com](mailto:info@adldentalab.com)



[adldentalab.com](http://adldentalab.com)

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China