



China Dental Lab Over 13 Years In Dental Lab Industry Providing Emax Veneer 16 Shades To Meet Diverse Dental Restoration Needs

Our Product Introduction

for more products please visit us on addentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians




Product Specification

- Experience: Over 13 Years In Dental Lab Industry
- Pricing model: Competitive Pricing With Volume Discounts
- Way Of Manufacture: Press And Porcelain Layered
- Esthetics: Natural
- Materials Used: BPS Acrylic/VITA Teeth
- Turnaround Time: In 72 Hours
- Product Type: Emax Veneer
- Service Type: China Dental Laboratory Services
- Highlight: **16 Shades eMax Veneer ,
Dental Lab eMax Veneer ,
Dental Restoration eMax Veneer**



More Images



Product Description

Our F

Product Description:

The China Dental Lab-Advanced Dental Lab offers exceptional dental prosthetics crafted with precision and care, making it an ideal choice for dental professionals seeking high-quality solutions. Utilizing premium materials such as BPS Acrylic and VITA Teeth, this product ensures durability, esthetic appeal, and patient comfort. The combination of these materials delivers a natural appearance and long-lasting performance, which is essential for successful dental restorations.

Manufactured using advanced techniques that include the Press and Porcelain Layered methods, the China Dental Lab product stands out for its superior craftsmanship. The Press technique ensures a strong, dense base for the dental prosthetics, while the porcelain layering adds a lifelike translucency and shade depth that closely mimics natural teeth. This meticulous manufacturing process results in prosthetics that not only fit perfectly but also enhance the overall smile aesthetics.

One of the standout features of this China Dental Laboratory Service is its fast turnaround time. With a commitment to efficiency and quality, the lab guarantees delivery within 72 hours, making it an excellent option for dental practices that require quick and reliable dental solutions. This fast turnaround time allows dental professionals to meet patient needs promptly without compromising on the quality of the final product.

The product is available in 16 shades, offering a broad spectrum of color options to match various natural tooth colors. This extensive shade selection ensures that dental restorations blend seamlessly with the patient's existing dentition, providing a natural and harmonious appearance. Whether the patient requires brighter shades or more subdued tones, the China Dental Lab can accommodate these preferences with precision.

Outsourcing to this China Dental Lab provides dental practices with access to top-tier dental laboratory services without the need for in-house manufacturing. By choosing to outsource dental lab work, dental professionals can benefit from specialized expertise, advanced manufacturing technologies, and cost-effective solutions. This service model not only saves time but also reduces overhead costs, allowing practices to focus more on patient care.

In summary, the China Dental Lab product offers a comprehensive solution for dental restorations, combining premium materials, sophisticated manufacturing techniques, and a wide range of shade options. With its fast turnaround time and reliable service, it is an excellent choice for dental professionals looking to outsource dental lab work. The use of BPS Acrylic and VITA Teeth materials, along with the Press and Porcelain Layered manufacturing methods, ensures that every prosthetic meets high standards of quality and esthetics. This product exemplifies the best in China Dental Laboratory Services, delivering both efficiency and excellence to dental practitioners worldwide.

Features:

Product Name: China Dental Lab

Materials Used: BPS Acrylic, VITA Teeth

Esthetics: Natural appearance for a flawless smile

Way of Manufacture: Press and porcelain layered technique

Shipping Options: UPS, DHL, FEDEX for reliable delivery

Service Type: China Dental Laboratory Services

Specialized in all on 6 dental solutions

Recognized as a top 5 dental lab for quality and precision

Expertise in all on 6 implant restorations

Technical Parameters:

Esthetics	Natural
Way Of Manufacture	Press And Porcelain Layered
Shipping options	UPS/DHL/FEDEX
Shade	16 Shades
Service type	China Dental Laboratory Services
Materials used	BPS Acrylic/VITA Teeth
Experience	Over 13 Years In Dental Lab Industry
Pricing model	Competitive Pricing With Volume Discounts
Turnaround Time	In 72 Hours
Product type	Emax Veneer

Applications:

ADL, a leading brand in the dental industry originating from China, stands out as one of the top 5 dental lab providers known for its exceptional quality and reliability. With over 13 years of experience in the dental lab industry, ADL specializes in delivering superior dental prosthetics using high-quality materials such as BPS Acrylic and VITA Teeth. These materials ensure a natural esthetic, making the final restorations blend seamlessly with the patient's natural dentition.

The ADL China Dental Lab product is ideal for a wide range of application occasions and scenarios. It is particularly suited for complex restorative cases such as the popular all on 6 dental implants, where precision and durability are paramount. Thanks to ADL's advanced manufacturing process, which combines press and porcelain layered techniques, dental professionals can trust that each prosthetic will meet the highest standards of craftsmanship and function.

One of the key strengths of ADL is its fast turnaround time, with delivery typically completed within 72 hours. This rapid delivery ensures that dental clinics can maintain efficient workflows and provide timely treatments to their patients. The supply ability of ADL is robust, supported by a skilled team of 630 technicians who are dedicated to producing top-quality dental restorations. This capacity allows for flexible minimum order quantities, which are negotiable to accommodate the specific needs of each client.

ADL's products are certified by ISO, CE, and FDA, guaranteeing compliance with international quality and safety standards. Packaging is thoughtfully designed, using plastic bags and boxes to ensure the products are protected during transit. Shipping options include reliable carriers such as UPS, DHL, and FEDEX, providing secure and timely delivery worldwide.

Payment terms are flexible, accepting Paypal or T/T, making transactions convenient for dental professionals globally. Whether you are a dental clinic seeking dependable all on 6 solutions or looking for a dental lab that offers a fast turnaround time without compromising quality, ADL China Dental Lab is an excellent choice. Its combination of advanced manufacturing, premium materials, certified quality, and efficient logistics makes it a trusted partner for dental professionals around the world.

Packing and Shipping:

Our China Dental Lab products are carefully packaged to ensure maximum protection during transit. Each item is securely placed in custom-fitted packaging materials that prevent movement and damage. The packaging is designed to maintain the integrity and hygiene of the dental products, using sterile, high-quality materials.

For shipping, we offer reliable and timely delivery services worldwide. Products are shipped via trusted carriers with full tracking capabilities. We also provide options for expedited shipping to meet urgent needs. All shipments comply with international shipping regulations and include necessary documentation to ensure smooth customs clearance.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China