



China Dental Lab Full Zirconia dental restorations dental lab specializing in AIDITE ZIRCONIA BLOCK materials for durable dental solutions

Our Product Introduction

for more products please visit us on adldentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Service Type: China Dental Laboratory Services
- Way Of Manufacture: Milled
- Shipping Options: UPS/DHL/FEDEX
- Esthetics: Natural
- Pricingmodel: Competitive Pricing With Volume Discounts
- Product Type: Full Zirconia
- Materials Used: AIDITE ZIRCONIA BLOCK
- Experience: Over 13 Years In Dental Lab Industry
- Strength: High
- Highlight: dental restorations ZIRCONIA BLOCK, China Dental Lab Full Zirconia, durable ZIRCONIA BLOCK



More Images



Product Description

Our F

Product Description:

With over 13 years of experience in the dental lab industry, our China dental lab ---Advanced Dental Lab has established itself as a trusted partner for dental professionals worldwide. We pride ourselves on delivering high-quality dental restorations using advanced technology and premium materials, ensuring precision and durability in every case. Our extensive expertise enables us to cater to a wide range of dental needs, making us the ideal choice for those looking to outsource dental lab services with confidence.

One of the core strengths of our China dental lab is our commitment to using only the best materials available. We utilize AIDITE Zirconia Block, renowned for its exceptional strength, esthetics, and biocompatibility. This material is perfect for fabricating crowns, bridges, and implant restorations that require both robustness and natural appearance. By integrating AIDITE Zirconia into our workflow, we ensure that every restoration meets the highest standards of quality and patient satisfaction.

Our product offerings include a comprehensive selection of shades to match diverse patient requirements. With 16 shades available, dental professionals can achieve highly accurate color matching, resulting in restorations that blend seamlessly with natural teeth. This wide shade range supports various clinical cases, including complex all on 4 implant procedures, where esthetics and functionality are paramount.

Understanding the time-sensitive nature of dental treatments, our China dental lab guarantees a fast turnaround time, delivering completed restorations within 72 hours. This rapid processing time allows dentists and clinics to provide efficient care while minimizing patient wait times. Whether you require single-unit crowns or full-arch restorations like all on 4 solutions, our lab's streamlined processes and skilled technicians ensure timely delivery without compromising quality.

We offer flexible and reliable shipping options to accommodate the needs of our global clients. Our partnerships with leading courier services such as UPS, DHL, and FEDEX enable us to provide secure and prompt delivery to any destination. This logistical capability ensures that outsourcing dental lab work to our China dental lab is both convenient and dependable, supporting smooth clinical workflows across different regions.

Outsourcing dental lab services to our China dental lab means gaining access to cutting-edge dental technology, experienced professionals, and superior materials, all while optimizing cost efficiency. Our extensive background in the industry and dedication to quality make us a preferred choice for dental practices seeking to enhance their prosthetic offerings. Whether your focus is on traditional crowns, bridges, or advanced procedures like all on 4 implants, our lab delivers consistent results that meet clinical expectations.

By choosing our China dental lab for your outsourcing needs, you benefit from a partner who understands the intricacies of dental restoration fabrication. We prioritize communication, accuracy, and customer satisfaction, ensuring that every case is handled with the utmost care. Our commitment to excellence and fast turnaround times enable dental professionals to provide superior care and improve patient outcomes effectively.

In summary, our China dental lab combines over 13 years of industry experience, the use of premium AIDITE Zirconia Block materials, a wide range of 16 shades, and a guaranteed 72-hour turnaround time to deliver top-quality dental restorations. With reliable shipping options through UPS, DHL, and FEDEX, we offer a seamless outsourcing solution for dental labs worldwide. Whether you require standard restorations or complex all on 4 implant solutions, our lab is dedicated to supporting your practice with precision, speed, and superior craftsmanship.

Features:

Product Name: China Dental Lab

Product Type: Full Zirconia

Materials Used: AIDITE ZIRCONIA BLOCK

Esthetics: Natural appearance for a lifelike smile

Shipping Options: UPS, DHL, FEDEX for reliable delivery

Pricing Model: Competitive pricing with volume discounts available

Specialized in china dental lab services including all on x and all on 4 dental solutions

Technical Parameters:

Materials Used	AIDITE ZIRCONIA BLOCK
Turnaround Time	In 72 Hours
Esthetics	Natural
Shipping Options	UPS/DHL/FEDEX
Experience	Over 13 Years In Dental Lab Industry
Shade	16 Shades
Service Type	China Dental Laboratory Services
Way Of Manufacture	Milled
Product Type	Full Zirconia
Pricing Model	Competitive Pricing With Volume Discounts

Applications:

The ADL China Dental Lab product is perfectly designed for a wide range of dental application occasions and scenarios, making it an ideal choice for dental professionals seeking reliable and high-quality dental laboratory services. As a reputable brand originating from China, ADL has over 13 years of experience in the dental lab industry, offering expertise and precision in every case handled. The lab specializes in outsourcing dental lab solutions, providing dental clinics and practitioners with a seamless way to enhance their treatment offerings without compromising on quality or speed.

One of the most prominent application scenarios for the ADL product is the "All on 6" dental implant procedure. This complex treatment requires exceptional accuracy and durable materials, and ADL meets these demands by utilizing AIDITE Zirconia Block materials, known for their strength and biocompatibility. Dental professionals can confidently outsource dental lab work related to All on 6 implants, knowing that the products will meet stringent ISO, CE, and FDA certifications, ensuring patient safety and satisfaction.

The ADL China Dental Lab also excels in delivering a fast turnaround time, with a guaranteed delivery within 72 hours. This rapid processing capability is crucial for dental practices looking to provide timely patient care and reduce waiting times. With a supply ability of 630 skilled technicians, ADL handles a large volume of orders efficiently without sacrificing quality. This allows dental clinics to maintain smooth workflows and meet urgent cases with ease.

Pricing is another advantageous aspect, as ADL offers competitive pricing with volume discounts, making it an economical choice for dental labs and clinics looking to outsource dental lab services. The minimum order quantity is negotiable, providing flexibility for both small and large-scale dental practices. Packaging is carefully done using plastic bags and boxes to ensure safe delivery, and payment terms are convenient, accepting PayPal or T/T transactions.

In summary, the ADL China Dental Lab product is suitable for various dental application scenarios, from routine restorations to complex implant cases like All on 6. Its combination of fast turnaround time, high-quality certified materials, experienced technicians, and flexible pricing models make it a trusted partner for dental professionals worldwide who want to outsource dental lab work efficiently and effectively.



Advanced Dental Laboratory



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China