



## China Dental Lab Providing AIDITE ZIRCONIA BLOCK Materials Alongside Turnaround Time In 72 Hours For Dental Manufacturing

Our Product Introduction

for more products please visit us on [addentalab.com](http://addentalab.com)

### Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Way Of Manufacture: Milled
- Materials Used: AIDITE ZIRCONIA BLOCK
- Turnaround Time: In 72 Hours
- Product Type: Full Zirconia
- Esthetics: Natural
- Shade: 16 Shades
- Service Type: China Dental Laboratory Services
- Pricingmodel: Competitive Pricing With Volume Discounts
- Strength: High
- Highlight: **AIDITE ZIRCONIA BLOCK Materials, Dental Manufacturing ZIRCONIA BLOCK**



### More Images



### Product Description

Our F

## Product Description:

China Dental Lab -Advanced Dental Lab is a premier provider of high-quality dental laboratory services, recognized as one of the top 5 dental lab facilities in the industry. With over 13 years of extensive experience, this dental lab has built a solid reputation for delivering exceptional dental solutions, specializing in full zirconia restorations. Their commitment to excellence, advanced technology, and skilled craftsmanship make them an ideal partner for dental professionals seeking reliable and efficient dental lab services. One of the standout features of China Dental Lab is its competitive pricing model designed to accommodate dental practices of all sizes. By offering volume discounts, they ensure that clients receive the best value without compromising quality. This pricing strategy allows dental clinics and practitioners to optimize their budgets while accessing top-tier dental lab products, making China Dental Lab a preferred choice among dental professionals worldwide.

China Dental Lab's specialization in full zirconia restorations reflects their dedication to providing durable, aesthetically pleasing, and biocompatible dental solutions. Full zirconia crowns and bridges manufactured by the lab are known for their strength, longevity, and natural appearance, which significantly enhance patient satisfaction. As a trusted China dental laboratory, they employ cutting-edge techniques and materials to produce restorations that meet the highest industry standards.

The lab also offers comprehensive China dental laboratory services tailored to meet the diverse needs of dental practices. Whether it's single crowns, bridges, implants, or complex cases such as the "all on 6" treatment concept, China Dental Lab delivers precision-crafted restorations with remarkable accuracy. Their expertise in handling "all on 6" cases, which involve full-arch implant-supported prostheses, positions them among the top 5 dental lab providers specializing in implant solutions.

Turnaround time is a critical factor in dental lab services, and China Dental Lab excels in this area by guaranteeing delivery within 72 hours. This rapid turnaround ensures that dental practices can provide timely treatment to their patients without unnecessary delays, thereby improving workflow efficiency and patient experience. The lab's streamlined production process and dedicated team work diligently to meet tight deadlines while maintaining exceptional quality standards.

With over 13 years in the dental lab industry, China Dental Lab combines experience, innovation, and customer-focused service to support dental professionals worldwide. Their extensive industry knowledge allows them to understand the unique challenges faced by dental practitioners and provide customized solutions that enhance clinical outcomes. As one of the top 5 dental lab facilities, they continue to expand their capabilities and product offerings to meet evolving market demands.

In summary, China Dental Lab stands out as a reliable and efficient partner for dental professionals seeking high-quality full zirconia restorations and comprehensive dental laboratory services. Their competitive pricing with volume discounts, fast 72-hour turnaround time, and specialization in complex procedures like the all on 6 implant system underscore their commitment to excellence. Whether you are a dental clinic looking to enhance your restorative options or a dental professional aiming for top-notch lab support, China Dental Lab offers unmatched expertise and service quality as one of the top 5 dental lab choices in the market.

## Features:

Product Name: China Dental Lab  
Shipping Options: UPS, DHL, FEDEX  
Shade Options: 16 Shades available  
Experience: Over 13 Years in Dental Lab Industry  
Service Type: China Dental Laboratory Services  
Manufacturing Method: Milled for precision and quality  
Specialized in outsource dental lab solutions  
Supports all on 6 dental implant cases  
Reliable outsource dental lab with extensive industry expertise

## Technical Parameters:

Turnaround Time	In 72 Hours
Shipping Options	UPS/DHL/FEDEX
Product Type	Full Zirconia
Way Of Manufacture	Milled
Service Type	China Dental Laboratory Services
Experience	Over 13 Years In Dental Lab Industry
Materials Used	AIDITE ZIRCONIA BLOCK
Esthetics	Natural
Pricing Model	Competitive Pricing With Volume Discounts
Shade	16 Shades

## Applications:

The ADL dental lab, based in China, is a premier outsource dental lab known for its exceptional quality and reliability. Specializing in advanced dental prosthetics, ADL caters to a wide range of dental restoration needs, including complex cases such as all on 4 and all on 6 implant solutions. These procedures require precision and esthetic excellence, both of which ADL guarantees through its meticulous milled manufacturing process.

With ISO, CE, and FDA certifications, ADL ensures that every product meets the highest international standards. The use of premium AIDITE zirconia blocks allows the lab to produce restorations that are not only durable but also natural-looking, offering patients the best combination of strength and esthetics. Whether you are a dental professional looking to outsource dental lab work or a clinic seeking reliable partners, ADL's supply ability of 630 skilled technicians ensures prompt and efficient service.

ADL understands the importance of quick turnaround times in dental treatments. With a standard delivery and turnaround time of 72 hours, dentists can expect timely delivery of all on 4 or all on 6 prosthetics, helping streamline patient workflows and reduce waiting periods. The packaging is carefully done in plastic bags and boxes to ensure safe transport, with multiple shipping options available, including UPS, DHL, and FEDEX, catering to global clients.

The lab offers flexible payment terms, accepting Paypal or T/T, and negotiable minimum order quantities and pricing, making it accessible to a wide range of dental practices and labs. This flexibility, combined with advanced milling techniques and high-quality zirconia materials, positions ADL as a trusted choice for outsourcing dental lab work.

In summary, ADL's China-based outsource dental lab is ideal for dental professionals focusing on implant-supported restorations such as all on 4 and all on 6. Its commitment to quality, fast delivery, certified processes, and natural esthetic outcomes make it a valuable partner in modern dental care.



**Advanced Dental Laboratory**



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong  
Province, P.R.China