



China Dental Lab outsource dental lab Ti reinforcement bar with full zirconia providing competitive pricing with volume discounts tailored to dental professionals requiring and manufacturing

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Pricingmodel: Competitive Pricing With Volume Discounts
- Service Type: China Dental Laboratory Services
- Way Of Manufacture: Milled
- Product Type: Full Zirconia With Reinforced Ibar
- Shipping Options: UPS/DHL/FEDEX
- Materials Used: AIDITE ZIRCONIA BLOCK
- Shade: 16 Shades
- Experience: Over 13 Years In Dental Lab Industry
- System: Ibar



More Images



Product Description

Our F

Product Description:

With over 13 years of expertise in the dental lab industry, our China dental lab-Advanced Dental Lab stands as a trusted partner for dental professionals worldwide seeking high-quality, reliable solutions. Our extensive experience has enabled us to refine our processes and deliver exceptional dental restorations that meet the highest standards of precision and esthetics. By choosing to outsource dental lab services to us, you benefit from our deep industry knowledge, advanced technology, and a commitment to excellence that ensures every product we deliver enhances patient satisfaction.

One of the key advantages of working with our -Advanced Dental Lab stands is our fast turnaround time. We understand the critical importance of timely delivery in dental practices, which is why we guarantee a turnaround time of just 72 hours. This rapid processing capability allows dental clinics and professionals to maintain seamless workflows, reduce patient wait times, and improve overall treatment outcomes. Whether you need crowns, bridges, dentures, or implants, our efficient production system ensures your orders are completed and shipped promptly without compromising quality.

Esthetics is a fundamental aspect of dental restorations, and our lab excels in creating natural-looking dental prosthetics that blend seamlessly with patients' existing teeth. Utilizing the latest materials and advanced fabrication techniques, we produce restorations that mimic the color, texture, and translucency of natural teeth. Our skilled technicians pay meticulous attention to detail, ensuring each piece not only fits perfectly but also enhances the patient's smile with a realistic appearance. When you outsource dental lab work to our -Advanced Dental Lab stands, you can be confident that the esthetic results will meet or exceed your expectations.

We also understand that pricing is a crucial consideration for dental practices. Our China dental lab-Advanced Dental Lab stands offers a competitive pricing model designed to provide exceptional value without sacrificing quality. In addition to affordable base rates, we provide volume discounts to support growing practices and larger orders. This pricing strategy enables dental professionals to optimize their budgets while accessing premium dental restorations. By partnering with us, you gain access to cost-effective solutions that enhance your practice's profitability and patient care capabilities.

Shipping logistics are another important factor when outsourcing dental lab services, and we offer flexible and reliable shipping options to meet your needs. Our China dental lab-Advanced Dental Lab stands works with leading carriers such as UPS, DHL, and FedEx to ensure safe, secure, and timely delivery of your dental restorations. Whether you require expedited shipping or standard delivery, we tailor our shipping methods to fit your schedule and location. This comprehensive shipping network supports our promise of a 72-hour turnaround time, allowing you to receive your orders quickly and start patient treatments without delay.

Choosing to outsource dental lab work to our China dental lab Advanced Dental Lab stands means partnering with a provider committed to quality, speed, and customer satisfaction. Our blend of over 13 years of industry experience, fast turnaround time, natural esthetics, competitive pricing, and dependable shipping solutions makes us an ideal choice for dental professionals looking to elevate their practice. We invite you to experience the benefits of working with a leading dental lab that understands your needs and delivers outstanding results every time.

Features:

Product Name: Ti reinforcement bar with full zirconia

Product Type: Full Zirconia and ibar

Materials Used: AIDITE Zirconia Block

Shade Options: 16 Shades available

Turnaround Time: In 72 Hours

Pricing Model: Competitive Pricing with Volume Discounts

Recognized as one of the top 5 dental labs in China

Specialized in all on 4 dental solutions

High-quality China Dental Lab trusted by professionals worldwide

Technical Parameters:

Way Of Manufacture	Milled
Esthetics	Natural
Experience	Over 13 Years In Dental Lab Industry
Product Type	Full Zirconia and ibar
Shipping Options	UPS/DHL/FEDEX
Materials Used	AIDITE ZIRCONIA BLOCK
Shade	16 Shades
Service Type	China Dental Laboratory Services
Pricing Model	Competitive Pricing With Volume Discounts
Turnaround Time	In 72 Hours

This China Dental Lab product is ideal for all on x procedures, including all on 6 solutions, providing high-quality, natural esthetics and reliable performance in all on 6 cases.

Applications:

The ADL Ti reinforcement bar with full zirconia dental lab product, originating from China and certified by ISO, CE, and FDA, is designed to meet the diverse needs of dental professionals worldwide. This product is ideal for various dental restoration applications, particularly in scenarios requiring high precision and durability, such as all on 6 implant procedures. The use of AIDITE Zirconia Block as the primary material ensures superior strength and esthetic results, making it a preferred choice for full arch restorations and other complex dental prosthetics.

Dental clinics and dental lab professionals looking to outsource dental lab work will find ADL's offering highly advantageous due to its fast turnaround time of just 72 hours. This quick delivery capability supports urgent cases and enhances workflow efficiency, enabling dentists to provide timely treatment to their patients. The product is milled with advanced technology by a team of 630 skilled technicians, ensuring consistent quality and precision in every piece fabricated.

ADL's dental lab services are suitable for a wide range of applications including single crowns, bridges, and full arch restorations, especially in implant-supported cases like all on 6. The flexible minimum order quantity, which is negotiable, allows both small and large dental practices to benefit from competitive pricing models with volume discounts. This makes it easier for clinics to manage costs without

compromising on quality.

Packaging is thoughtfully handled with protective plastic bags and boxes to safeguard the zirconia restorations during shipping. Payment terms are flexible, accepting Paypal or T/T, which facilitates smooth international transactions. With a supply ability supported by a large technician workforce, ADL ensures reliable availability and consistent supply for dental labs looking to outsource their manufacturing needs.

In summary, ADL's Ti reinforcement bar with full zirconia product is perfect for dental professionals seeking a trustworthy outsource dental lab solution that offers competitive pricing, high-quality milled zirconia restorations, and a fast turnaround time. Whether for all on 6 implant cases or other restorative dental applications, ADL provides a reliable, certified, and efficient service tailored to the dynamic needs of modern dental practices.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China