



## Competitive Pricing With Volume Discounts Dental Lab Offering 16 Shades Dental Prosthetics Tailored To Medical Industry

Our Product Introduction

for more products please visit us on [adldentalab.com](http://adldentalab.com)

### Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Materials Used: AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
- Shipping Options: UPS/DHL/FEDEX
- Pricingmodel: Competitive Pricing With Volume Discounts
- Experience: Over 13 Years In Dental Lab Industry
- Shade: 16 Shades
- Service Type: China Dental Laboratory Services
- Inspection: High Magnification Inspection
- Way Of Manufacture: Milled

### Product Description

**Product Description:**

Our Product

When it comes to precision, reliability, and exceptional craftsmanship in dental restorations, the China Dental Lab stands out as a premier choice for dental professionals worldwide. With over 13 years of experience in the dental lab industry, this lab has consistently demonstrated its commitment to quality and innovation, making it one of the top 5 dental labs highly recommended by dentists and clinicians alike.

Specializing in Porcelain Fused To Zirconia (PFZ) restorations, the China Dental Lab offers products that combine the strength and durability of zirconia with the natural aesthetics of porcelain. This hybrid approach ensures that patients receive restorations that are not only robust and long-lasting but also visually appealing, matching the natural teeth in color and translucency. Whether you are looking for crowns, bridges, or full-arch restorations, including complex all on 4 solutions, this lab delivers exceptional results tailored to meet the unique needs of each patient.

One of the key advantages of partnering with the China Dental Lab is their seamless and efficient manufacturing process. Utilizing advanced milling technology, all products are milled with precision to ensure an exact fit and superior finish. This method allows for consistent quality and repeatability, which is crucial for dental professionals who require dependable restorations every time. The use of milled porcelain fused to zirconia also enhances the structural integrity of the restorations, reducing the risk of fractures and ensuring long-term success in clinical applications.

Efficiency is another hallmark of the China Dental Lab. Understanding the fast-paced nature of dental practices, the lab guarantees a turnaround time of just 72 hours. This rapid delivery ensures that dental clinics can maintain their schedules without compromising on quality, enabling quicker patient treatment and satisfaction. For dental professionals looking to outsource dental lab work, this quick turnaround is invaluable, helping to streamline operations and improve overall workflow.

Customer service is a cornerstone of the China Dental Lab's business philosophy. Providing support 7 days a week, the lab ensures that clients have access to assistance whenever needed. Whether you have questions about your order, require technical advice, or need help with case submissions, their dedicated customer service team is always available to provide prompt and professional support. This 7-day availability demonstrates their commitment to client satisfaction and ensures smooth communication throughout every step of the process.

Choosing to outsource dental lab work to the China Dental Lab means partnering with a trusted leader in the dental industry. Their extensive experience, combined with cutting-edge technology and a strong customer-centric approach, positions them as a top choice for dental professionals seeking high-quality porcelain fused to zirconia restorations. This lab's expertise in handling complex cases, including all on 4 implant-supported prosthetics, further underscores their capability to meet diverse clinical requirements.

In summary, the China Dental Lab offers a comprehensive solution for dental labs and clinics looking to outsource dental lab fabrication with confidence. With over 13 years of experience, advanced milled manufacturing processes, a rapid 72-hour turnaround, and exceptional customer service available 7 days a week, they have earned their place among the top 5 dental labs in the industry. Trust the China Dental Lab to deliver durable, aesthetically pleasing, and reliable porcelain fused to zirconia restorations that enhance patient outcomes and support the success of your dental practice.

## Features:

Product Name: China Dental Lab

Experience: Over 13 Years In Dental Lab Industry

Inspection: High Magnification Inspection for superior quality

Way Of Manufacture: Milled for precision and durability

Strength: High strength materials ensuring long-lasting results

Shade Options: Available in 16 Shades to match natural teeth

Recognized as one of the top 5 dental labs in the industry

Fast turnaround time to meet urgent dental restoration needs

Trusted China dental lab with a proven track record

## Technical Parameters:

Materials Used	AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
Shipping Options	UPS/DHL/FEDEX
Service Type	China Dental Laboratory Services
Strength	High
Target Customers	International Dental Lab And Dentist
Pricing Model	Competitive Pricing With Volume Discounts
Customer Service	7 days A Week
Turnaround Time	In 72 Hours
Experience	Over 13 Years In Dental Lab Industry
Shade	16 Shades

## Applications:

The ADL Porcelain Fused To Zirconia product, originating from China, is an ideal choice for international dental labs and dentists seeking high-quality dental restorations. As one of the top 5 dental labs globally, ADL prides itself on delivering exceptional craftsmanship with a fast turnaround time of just 72 hours. This rapid delivery ensures that dental professionals can meet their patients' needs efficiently without compromising on quality.

Designed with versatility in mind, ADL's porcelain fused to zirconia restorations are perfectly suited for a wide range of dental applications, including complex procedures such as all on x treatments. The product is milled with precision by a team of 630 skilled technicians, ensuring consistency and durability in every piece. With 16 available shades, dental professionals can easily match the restorations to the natural teeth, providing an aesthetically pleasing and functional outcome for patients.

ADL adheres to stringent quality standards, holding ISO, CE, and FDA certifications, which guarantees that every product meets international safety and performance requirements. The manufacturing process is carefully controlled to maintain the highest quality, and the product is securely packaged in plastic bags and boxes to protect it during shipping.

Flexibility is a key attribute of ADL's service. The minimum order quantity is negotiable, allowing dental labs and dentists to order according to their specific needs. Pricing is also negotiable, making it accessible for a wide range of customers. Payment can be made conveniently via Paypal or T/T, facilitating smooth transactions across borders.

Whether you are working on individual crowns, bridges, or full-arch restorations like all on x, ADL's porcelain fused to zirconia products provide reliability, precision, and an excellent aesthetic outcome. The combination of advanced milling technology, expert technicians, and a commitment to fast turnaround time makes ADL a trusted partner for dental professionals worldwide.



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong  
Province, P.R.China