



China Dental Lab Porcelain Fused To Zirconia High Strength Providing Tailored Services To International Dental Lab And Dentist

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Esthetics: Natural
- Turnaround Time: In 72 Hours
- Shipping Options: UPS/DHL/FEDEX
- Product Type: Porcelain Fused To Zirconia
- Strength: High
- Inspection: High Magnification Inspection
- Way Of Manufacture: Milled
- Customer Service: 7days A Week



More Images



Product Description

Product Description:

Our F

China Dental Lab--Advanced Dental Lab is a premier provider of high-quality dental products and services, boasting over 13 years of extensive experience in the dental lab industry. Renowned as one of the top 5 dental labs globally, ADL has built a strong reputation for delivering exceptional dental restorations that meet the precise needs of international dental labs and dentists. With a focus on innovation, accuracy, and reliability, this lab specializes in producing milled dental prosthetics that ensure superior fit and durability. One of the standout features of ADL is its dedication to catering to the evolving demands of dental professionals worldwide. The lab offers a comprehensive range of dental solutions, including the highly sought-after all on x treatment options. This advanced approach to full-arch restoration has gained immense popularity among dental practitioners, and ADL's expertise in fabricating milled restorations for all on x cases ensures optimal patient outcomes.

ADL provides an extensive selection of 16 shades, allowing dentists to achieve the most natural and aesthetically pleasing results for their patients. This wide shade range facilitates perfect color matching, which is crucial for anterior and posterior restorations alike. Whether working with single crowns, bridges, or complex full-arch prosthetics, clients can rely on the lab's precision milling techniques and superior materials to deliver restorations that blend seamlessly with the patient's natural dentition.

Understanding the importance of timely delivery in the dental industry, China Dental Lab-ADL offers flexible and reliable shipping options via UPS, DHL, and FEDEX. This ensures that international dental labs and dentists receive their orders promptly and securely, no matter where they are located. The lab's streamlined logistics and global shipping network contribute to its status as a trusted partner for dental professionals across continents.

China Dental Lab-ADL's manufacturing process centers on advanced milling technology, which enhances the accuracy and consistency of each dental restoration. Milled restorations from China Dental Lab are known for their precision fit, strength, and esthetic appeal. By utilizing computer-aided design and computer-aided manufacturing (CAD/CAM) systems, the lab minimizes human error and maximizes efficiency, enabling the production of highly customized prosthetics tailored to individual patient needs.

For international dental labs and dentists seeking a reliable partner capable of delivering high-quality milled restorations with a broad shade range and excellent service, China Dental Lab-ADL stands out as a top choice. Their commitment to quality, innovation, and customer satisfaction has earned them a distinguished place among the top 5 dental labs worldwide. With over 13 years of experience, China Dental Lab continues to push the boundaries of dental technology and craftsmanship, making it an indispensable resource for dental professionals aiming to provide outstanding care to their patients.

In summary, China Dental Lab-ADL combines decades of industry expertise, advanced milling manufacturing methods, a diverse palette of 16 shades, and comprehensive shipping options including UPS, DHL, and FEDEX, to serve the needs of international dental labs and dentists. Their specialization in all on x restorations further highlights their capability to handle complex dental cases with precision and excellence. Choose China Dental Lab for dependable, high-quality dental prosthetics crafted to the highest standards in the dental industry.

Features:

Product Name: Porcelain Fused To Zirconia
 Product Type: Outsourcing service to international dentists and dental lab
 Shipping Options: UPS, DHL, FEDEX
 Turnaround Time: In 72 Hours
 Esthetics: Natural appearance for a seamless smile
 Strength: High durability for long-lasting use
 Recognized as a Top 5 dental lab for quality and service

Technical Parameters:

Service Type	China Dental Laboratory Services
Esthetics	Natural
Target Customers	International Dental Lab And Dentist
Strength	High
Experience	Over 13 Years In Dental Lab Industry
Shipping Options	UPS / DHL / FEDEX
Inspection	High Magnification Inspection
Shade	16 Shades
Customer Service	7 days A Week
Product Name	Porcelain Fused To Zirconia

Applications:

ADL, a renowned China Dental Lab brand, offers exceptional dental laboratory services that cater to a wide range of application occasions and scenarios in the dental industry. With its origin in China and certifications including ISO, CE, and FDA, ADL ensures high-quality standards suitable for dental professionals worldwide. The lab's fast turnaround time of 72 hours makes it an ideal partner for clinics and dental practices requiring quick and reliable solutions without compromising on quality.

One of the primary application scenarios for ADL's dental lab services is the fabrication of "All on 4" dental implant prosthetics. This revolutionary treatment requires precision and speed to restore patients' smiles efficiently. ADL's expert technicians, totaling 630 in number, utilize advanced technology and competitive pricing models with volume discounts to provide cost-effective solutions for large dental practices or individual dentists. The lab supports a minimum order quantity that is negotiable, allowing flexibility for both small and large-scale orders.

ADL's packaging details include secure plastic bags and boxes, ensuring that all dental products arrive safely and intact. With multiple shipping options such as UPS, DHL, and FedEx, the products can be delivered swiftly and reliably to clients around the globe. The combination of fast turnaround time and reliable shipping makes ADL a preferred choice for urgent dental cases requiring immediate attention.

The lab's customer service operates seven days a week, providing continuous support to dental professionals and ensuring smooth communication throughout the ordering and delivery process. Payment terms are flexible, accepting both PayPal and T/T, accommodating various client preferences. ADL's commitment to quality, speed, and excellent customer service positions it as a leading China dental laboratory service provider for diverse dental restoration needs.

Whether it is routine crown and bridge work, complex implant cases such as All on 4, or customized dentures, ADL's dental lab services cover a broad spectrum of dental restoration scenarios. Clients benefit from the reliable supply ability, competitive pricing, and the assurance of receiving products certified to meet international standards. This makes ADL an invaluable partner for dental professionals seeking efficiency, quality, and affordability in their dental lab collaborations.



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China