



Competitive Pricing With Volume Discounts outsource dental lab offering expert dental appliance fabrication with fast delivery options

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Shipping Options: UPS/DHL/FEDEX
- Turnaround Time: In 72 Hours
- Target Customers: International Dental Lab And Dentist
- Shade: 16 Shades
- Way Of Manufacture: Milled
- Pricingmodel: Competitive Pricing With Volume Discounts
- Materials Used: AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
- Strength: High

Product Description

Product Description:

Our Product

The China Dental Lab is a premier outsource dental lab renowned for its commitment to quality, precision, and customer satisfaction. Specializing in dental restorations crafted from premium materials, this lab uses the highly acclaimed AIDITE ZIRCONIA BLOCK to ensure that each product meets the highest standards of durability and esthetics. The AIDITE ZIRCONIA BLOCK is both CE and FDA approved, providing dental professionals and patients with the assurance that the materials used comply with stringent international safety and quality regulations. This focus on certified materials guarantees that every restoration produced is not only strong and reliable but also safe for long-term use.

One of the standout features of the China Dental Lab is its dedication to delivering natural esthetics in every dental piece. Whether you require crowns, bridges, or full arch restorations such as the popular all on x solutions, the lab's skilled technicians meticulously craft each product to mimic the natural appearance of teeth. The result is restorations that blend seamlessly with the patient's existing dentition, enhancing both confidence and comfort. The natural esthetics offered by the lab make it an ideal choice for dental professionals who want to provide their patients with restorations that look and feel like real teeth.

Strength is a critical factor in dental restorations, and the China Dental Lab excels in this area by utilizing the high-strength properties of the AIDITE ZIRCONIA BLOCK. This material is known for its exceptional toughness and resistance to fracture, making it suitable for a wide range of dental applications, including demanding full arch restorations like all on x treatments. Patients benefit from restorations that withstand the rigors of daily use, ensuring longevity and reducing the need for frequent repairs or replacements. The lab's focus on strength and durability reflects its commitment to producing reliable dental solutions that meet the needs of both dentists and patients. In addition to superior materials and craftsmanship, the China Dental Lab offers flexible and efficient shipping options to accommodate the needs of dental professionals worldwide. Clients can choose from trusted courier services such as UPS, DHL, and FEDEX to ensure timely and secure delivery of their dental restorations. This reliable shipping network helps streamline the workflow for dental practices, allowing them to outsource their dental lab needs with confidence that their orders will arrive promptly and in excellent condition.

Customer service is a cornerstone of the China Dental Lab's business philosophy. Understanding the importance of clear communication and support in the dental industry, the lab provides customer service seven days a week. This availability ensures that dental professionals can receive assistance, ask questions, and resolve any issues quickly and efficiently, minimizing downtime and enhancing overall satisfaction. The lab's commitment to excellent customer service reinforces its position as a trusted partner for dental practices looking to outsource dental lab work.

Choosing to outsource dental lab services to the China Dental Lab means partnering with a provider that values quality, reliability, and customer care. Whether you are seeking solutions for single-unit restorations or comprehensive all on x cases, the lab's expertise and use of CE and FDA approved materials guarantee superior outcomes. Their natural esthetic designs combined with the high strength of AIDITE ZIRCONIA BLOCK ensure that patients receive restorations that are both beautiful and durable. Flexible shipping options through UPS, DHL, and FEDEX, along with dedicated customer service seven days a week, make this lab an excellent choice for dental professionals worldwide.

In summary, the China Dental Lab stands out as a leading outsource dental lab that delivers high-quality zirconia restorations using CE and FDA approved AIDITE ZIRCONIA BLOCK material. The lab's focus on natural esthetics, exceptional strength, and comprehensive customer support ensures that dental professionals can confidently rely on their services. With convenient shipping options and a commitment to excellence, the China Dental Lab is your ideal partner for all your dental restoration needs, including complex all on x procedures.

Features:

Product Name: China Dental Lab

Target Customers: International Dental Lab and Dentist

Experience: Over 13 Years in Dental Lab Industry

Product Type: Porcelain Fused to Zirconia

Strength: High strength for durable dental restorations

Esthetics: Natural appearance for a seamless smile

Specialized in china dental lab services for all on 4 treatments

Fast turnaround time to meet urgent dental lab needs

Technical Parameters:

Customer Service	7 days A Week
Pricing Model	Competitive Pricing With Volume Discounts
Target Customers	International Dental Lab And Dentist
Way Of Manufacture	Milled
Experience	Over 13 Years In Dental Lab Industry
Strength	High
Turnaround Time	In 72 Hours
Esthetics	Natural
Materials Used	AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
Product Type	Porcelain Fused To Zirconia

Applications:

The ADL China dental lab is a leading provider in the dental industry, renowned for its porcelain fused to zirconia products. With over 13 years of experience in the dental lab industry, ADL has established itself as one of the top 5 dental labs, trusted by international dental labs and dentists worldwide. The brand's commitment to quality is reflected in its ISO, CE, and FDA certifications, ensuring that every product meets stringent global standards.

ADL specializes in a variety of dental restoration solutions, making it an ideal choice for numerous application occasions and scenarios. One of the most significant applications is for the "All on X" dental implant procedure. This innovative technique requires highly durable and aesthetic dental prosthetics, which ADL provides with its porcelain fused to zirconia products. These restorations offer excellent strength, biocompatibility, and natural appearance, making them perfect for full-arch implant restorations.

International dentists and dental labs rely on ADL for their complex cases, including crowns, bridges, and implant-supported restorations. The lab's ability to handle customized orders with a negotiable minimum order quantity makes it flexible for both small clinics and large dental laboratories. Moreover, the packaging details, consisting of secure plastic bags and boxes, ensure safe delivery of delicate dental prosthetics worldwide.

ADL's supply ability is supported by a skilled team of 630 technicians, enabling a fast turnaround time of just 72 hours. This quick delivery time is crucial for dental professionals who need reliable and timely solutions for their patients. Customers can place orders with confidence, knowing that payment terms are flexible, accepting both Paypal and T/T methods, and customer service is available 7 days a week to provide support throughout the ordering and delivery process.

In summary, the ADL China dental lab is perfectly suited for a wide range of dental restoration needs, especially in complex scenarios such as the all on x procedure. Its status as one of the top 5 dental labs globally, combined with its certifications, experienced technicians, and efficient logistics, makes it a preferred partner for international dental labs and dentists seeking high-quality porcelain fused to zirconia products.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China