



16 Shades Dental Lab Featuring Professional Dental Laboratory Services 7days A Week Precision and Reliability Guaranteed

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Strength: High
- Customer Service: 7days A Week
- Turnaround Time: In 72 Hours
- Way Of Manufacture: Milled
- Inspection: High Magnification Inspection
- Shade: 16 Shades
- Pricingmodel: Competitive Pricing With Volume Discounts
- Target Customers: International Dental Lab And Dentist

Product Description

Product Description:

Our Product

China Dental Lab offers a premium Porcelain Fused To Zirconia product designed specifically for international dental labs and dentists seeking high-quality dental restorations. Our product stands out for its exceptional esthetics, providing a natural appearance that closely mimics the translucency and shade of real teeth. This ensures that patients receive restorations that are not only durable but also visually pleasing, enhancing their smile with confidence.

Crafted using AIDITE ZIRCONIA BLOCK, a material renowned for its superior strength and biocompatibility, our porcelain fused to zirconia restorations guarantee long-lasting performance. The AIDITE zirconia is CE and FDA approved, ensuring that every product meets stringent international safety and quality standards. This compliance underscores our commitment to delivering products that dental professionals can trust for their patients' health and satisfaction.

Our China dental lab specializes in serving the needs of international dental labs and dentists, understanding the importance of reliability and precision in dental restorations. We cater to a broad range of dental restoration needs, including the increasingly popular all on 6 procedure. By using our porcelain fused to zirconia solutions, dental professionals can achieve optimal results for patients undergoing full-arch implant restorations, combining strength, esthetics, and functionality.

One of the key advantages of working with our China dental lab is our fast turnaround time. We recognize that timely delivery is crucial in dental practice to ensure patient satisfaction and efficient workflow. Our streamlined production process and expert craftsmanship allow us to deliver high-quality restorations quickly without compromising on quality, helping dental professionals meet tight deadlines and urgent case requirements.

Customer service is at the heart of our operations. We offer support 7 days a week to assist international clients with orders, technical inquiries, and after-sales service. Our dedicated customer service team is always ready to provide guidance, ensuring a smooth and pleasant experience from order placement to delivery. This commitment to client satisfaction reinforces our position as a reliable partner in the dental industry.

In summary, our porcelain fused to zirconia product from the China dental lab combines natural esthetics, durable materials, and compliance with international standards, making it an ideal choice for dental professionals worldwide. Whether you are performing routine restorations or complex all on 6 implant procedures, our products provide the quality and reliability you need. With fast turnaround times and comprehensive customer support available seven days a week, we ensure that your practice runs smoothly and your patients receive the best possible care.

Choose our China dental lab for your porcelain fused to zirconia needs and experience the perfect blend of advanced materials, expert craftsmanship, and dedicated service tailored to the global dental community. Elevate your dental restorations with products that deliver natural beauty, strength, and precision, all backed by our commitment to excellence and timely delivery.

Features:

Product Name: China Dental Lab

Inspection: High Magnification Inspection for precise quality control

Product Type: Porcelain Fused To Zirconia

Target Customers: International Dental Lab and Dentist

Strength: High strength for durable dental restorations

Materials Used: AIDITE Zirconia Block; CE and FDA Approved Material ensuring safety and reliability

Fast turnaround time, delivering efficient service to meet urgent needs

Reliable china dental lab known for quality and quick delivery

Technical Parameters:

Shipping options	UPS/DHL/FEDEX
Strength	High
Experience	Over 13 Years In Dental Lab Industry
Product type	Porcelain Fused To Zirconia
Shade	16 Shades
Materials used	AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
Pricing model	Competitive Pricing With Volume Discounts
Service type	China Dental Laboratory Services
Inspection	High Magnification Inspection
Target customers	International Dental Lab And Dentist

Applications:

The ADL China Dental Lab is a premier choice for dental professionals seeking high-quality, reliable dental prosthetics and restorations. With its advanced milled manufacturing process, ADL ensures each product boasts natural esthetics that meet the highest standards of dental craftsmanship. This makes it an excellent solution for a wide range of dental applications, including the popular all on 4 treatment concept, where precision and durability are paramount.

Dental clinics and laboratories rely on ADL China Dental Lab for its fast turnaround time, with a guaranteed delivery window of just 72 hours. This rapid delivery supports urgent cases and helps dental professionals maintain efficient workflows without compromising on quality. Whether you are handling routine restorations or complex full-arch rehabilitations like all on 4 implants, ADL's competitive pricing with volume discounts makes it a cost-effective partner, especially when minimum order quantities are negotiable to fit the needs of any sized practice.

Thanks to its ISO, CE, and FDA certifications, the ADL brand assures compliance with global quality and safety standards, making it suitable for a diverse range of clinical scenarios. Dental labs and clinics benefit from a supply ability supported by 630 skilled technicians, ensuring consistent quality and reliability for every order. The packaging is designed for safety and convenience, using plastic bags and boxes to protect delicate dental components during shipping.

ADL's customer service operates seven days a week, providing continuous support to address any inquiries or customized requirements. The lab offers flexible payment terms via Paypal or T/T and reliable shipping options through UPS, DHL, and FedEx, facilitating smooth international deliveries. This comprehensive service model makes ADL China Dental Lab an ideal partner for dental professionals focusing on efficiency, quality, and esthetic excellence in treatments such as all on 4 and beyond.

In summary, ADL China Dental Lab is perfectly suited for dental professionals who demand natural-looking restorations, quick turnaround times, and dependable service. Its advanced milled manufacturing, competitive pricing, and global certifications make it a standout option for all dental restoration needs, especially in fast-paced clinical environments where time and quality are critical.



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China