



## Top 5 China dental lab providers Porcelain Fused To Zirconia ensuring turnaround time in 72 hours service for dental offices and labs

Our Product Introduction

for more products please visit us on [adidentalab.com](http://adidentalab.com)

### Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Materials Used: AIDITE ZIRCONIA BLOCK; CE And FDA Approved Material
- Shipping Options: UPS/DHL/FEDEX
- Way Of Manufacture: Milled And Layered
- Turnaround Time: In 72 Hours
- Pricingmodel: Competitive Pricing With Volume Discounts
- Product Type: Porcelain Fused To Zirconia
- Target Customers: International Dental Lab And Dentist
- Esthetics: Natural



### More Images



### Product Description

**Product Description:**

Our F

Our China Dental Lab---Advanced Dental Lab is a premier outsource dental lab that specializes in providing high-quality dental restorations with a fast turnaround time. We understand the critical importance of efficiency and precision in dental procedures, which is why our lab is dedicated to delivering superior products within 72 hours. Whether you are working on single crowns, bridges, or complex All on 6 cases, our lab ensures you receive reliable and timely dental solutions that meet your clinical and aesthetic requirements. One of the key features that set our China Dental Lab apart is our commitment to using only the best materials available. We utilize AIDITE Zirconia Block, a state-of-the-art ceramic material known for its outstanding strength, durability, and biocompatibility. This material is CE and FDA approved, guaranteeing that all dental restorations produced in our lab meet the highest international standards for safety and quality. By selecting AIDITE Zirconia, we ensure that your dental prosthetics not only look natural but also provide long-lasting performance for your patients.

The manufacturing process in our lab is fully digital and employs advanced milling technology. Our skilled technicians use precision milling machines to craft each restoration with meticulous attention to detail. This method allows us to produce dental restorations with excellent fit and finish, minimizing adjustments and chair time for dental professionals. The milling process also contributes to the lab's ability to maintain a fast turnaround time, enabling you to receive your orders within 72 hours consistently.

Understanding the diverse needs of dental professionals, our China Dental Lab offers an extensive shade range with 16 shades available. This wide selection allows you to match the dental restoration perfectly to the patient's natural teeth, enhancing the overall esthetic outcome. Whether you are restoring anterior or posterior teeth, our shade options provide the flexibility required to achieve optimal results for every case.

Customer service is a cornerstone of our operation. We provide support 7 days a week to ensure seamless communication and assistance whenever you need it. Our dedicated customer service team is ready to help with order inquiries, technical questions, and custom requests, making your experience with our outsource dental lab smooth and hassle-free. This level of accessibility demonstrates our commitment to supporting dental professionals throughout their workflow.

Partnering with our China Dental Lab means entrusting your dental cases to a reliable and efficient provider who prioritizes quality, speed, and customer satisfaction. Whether you specialize in All on 6 implant restorations or require precision-milled crowns and bridges, our lab is equipped to deliver exceptional products that meet your clinical demands. Our fast turnaround time of 72 hours ensures that you can keep your practice running smoothly without compromising on quality.

In summary, our China Dental Lab offers an unbeatable combination of premium materials, advanced manufacturing methods, extensive shade options, and exemplary customer service. With AIDITE Zirconia Blocks that are CE and FDA approved, a fast turnaround time, and 7-day-a-week support, we provide dental professionals with the confidence and convenience they need to excel. Outsource your dental lab needs to us and experience the difference of a partner committed to your success and patient satisfaction.

## Features:

Product Name: Porcelain Fused To Zirconia

Esthetics: Natural appearance for a flawless smile

Customer Service: Available 7 days a week to assist you

Shipping Options: Reliable delivery via UPS, DHL, and FEDEX

Materials Used: AIDITE Zirconia Block; CE and FDA approved materials ensuring safety and quality

Inspection: High magnification inspection for precise and accurate dental restorations

Recognized as one of the top 5 dental labs in the industry

Specialized expertise in all on 4 dental implant solutions

## Technical Parameters:

Target Customers	International Dental Lab And Dentist
Product Type	Porcelain Fused To Zirconia
Experience	Over 13 Years In Dental Lab Industry
Service Type	China Dental Laboratory Services
Inspection	High Magnification Inspection
Turnaround Time	In 72 Hours
Shade	16 Shades
Pricing Model	Competitive Pricing With Volume Discounts
Strength	High
Esthetics	Natural

## Applications:

ADL, a renowned dental lab brand originating from China, stands out as one of the top 5 dental labs globally, trusted by dental professionals for its exceptional quality and reliability. With over 13 years of experience in the dental lab industry, ADL specializes in providing premium dental solutions, including the highly sought-after AIDITE ZIRCONIA BLOCK made from CE and FDA approved materials. These materials ensure high strength and natural esthetics, making ADL an ideal choice for a wide range of dental applications.

One of the main application occasions for ADL's products is in the provision of all on x dental restorations. This cutting-edge procedure requires precision, durability, and a natural look, all of which ADL's zirconia blocks deliver effectively. Dental practitioners who outsource dental lab work to ADL benefit from the lab's competitive pricing model, which includes volume discounts, making it cost-effective without compromising quality. This makes ADL a preferred partner for dental clinics seeking reliable, high-quality dental restorations within a quick turnaround time.

ADL's supply ability is supported by a skilled workforce of 630 technicians, ensuring consistent quality and timely delivery. The packaging details include secure plastic bags and boxes, maintaining the integrity of each product during shipment. With a delivery time of just 72 hours and flexible payment terms via Paypal or T/T, ADL offers an efficient and customer-friendly purchasing experience.

Whether dental professionals are looking to outsource dental lab services for routine crowns, bridges, or complex full-mouth restorations such as all on x, ADL provides a comprehensive solution that meets international certification standards including ISO, CE, and FDA. The negotiable minimum order quantity and pricing allow for tailored orders suited to the needs of different dental practices, from small clinics to large dental chains.

In summary, ADL from China represents a reliable and innovative dental lab partner. Its high-strength zirconia materials, natural esthetics, experienced technicians, and commitment to quality position it as a leading choice for dental professionals worldwide who want to outsource dental lab work with confidence and efficiency.





+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong  
Province, P.R.China