



## Long-lasting And Wear-resistant Crown And Bridge Dental Lab Offering Dental Laboratory Service Alongside 724H Customer Service

Our Product Introduction

for more products please visit us on [adidentalab.com](http://adidentalab.com)

### Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



### Product Specification

- Producttype: Dental Laboratory Service
- Hypoallergenic: Yes
- Experience : More Than 13 Years
- Shade: A1/A2/A3/A3.5...16 Shades
- Flexural Strength: 360 MPa
- Shipping: DHL/UPS/FEDEX
- Productname: Crown And Bridge Dental Lab
- Chemical Solubility: < 100 µg/mm²

### Product Description

**Product Description:**

Our Product

The Crown And Bridge Dental Lab product is a premium solution designed to meet the highest standards in restorative dentistry, offering exceptional quality and reliability for dental professionals and patients alike. This product is meticulously crafted to provide outstanding performance, combining advanced materials and expert craftsmanship to ensure optimal results in both function and aesthetics. One of the key features of the Crown And Bridge Dental Lab product is its use of porcelain fused to zirconia technology. This combination creates restorations that are not only strong and durable but also exhibit a natural-looking appearance that closely mimics real teeth. The porcelain layer enhances the esthetic appeal, providing a high esthetic finish that satisfies the most discerning patients and dental practitioners. The layered zirconia structure delivers an ideal balance between strength and beauty, ensuring that the crowns and bridges can withstand the rigors of daily use while maintaining their attractive appearance over time.

In addition to its superior esthetics, this dental lab product is designed with patient safety and comfort in mind. It is hypoallergenic and biocompatible, making it an excellent choice for individuals with sensitivities or allergies to certain dental materials. The biocompatibility of the materials used ensures that the restorations integrate seamlessly with the natural tissues in the mouth, minimizing the risk of adverse reactions and promoting long-term oral health. This focus on patient well-being enhances the overall treatment experience and contributes to the success of restorative procedures.

Strength and durability are paramount in any dental restoration, and the Crown And Bridge Dental Lab product excels in these areas. The high-strength zirconia framework provides exceptional resistance to fracture and wear, making it suitable for both anterior and posterior restorations. This durability ensures that the crowns and bridges maintain their structural integrity and function effectively over extended periods, reducing the need for replacements or repairs. Dental professionals can trust this product to deliver consistent, long-lasting results that meet the demands of everyday oral function.

Customer service is a cornerstone of the Crown And Bridge Dental Lab experience. To support dental practices worldwide, the company offers 7\*24H customer service, ensuring that assistance is available at any time of day or night. This round-the-clock support helps address any questions or concerns promptly, facilitating smooth communication and efficient order processing. Whether it's technical guidance, product inquiries, or after-sales support, the dedicated customer service team is committed to providing professional and timely responses.

Efficient and reliable shipping is another important aspect of this product offering. The Crown And Bridge Dental Lab product is shipped via trusted courier services such as DHL, UPS, and FEDEX. These global carriers guarantee fast and secure delivery, allowing dental professionals to receive their orders promptly and without compromise. Timely shipping helps maintain the workflow in dental clinics and labs, ensuring that restorative treatments can proceed without unnecessary delays.

In summary, the Crown And Bridge Dental Lab product stands out as an exceptional choice for dental restorations requiring high esthetic value, strength, and biocompatibility. Its porcelain fused to zirconia and layered zirconia construction offers a natural-looking finish that satisfies both aesthetic and functional demands. Coupled with hypoallergenic properties, 7\*24H customer service, and reliable shipping options through DHL, UPS, and FEDEX, this product provides a comprehensive solution for modern dental practices aiming to deliver superior care and patient satisfaction.

## Features:

Product Name: Crown And Bridge Dental Lab

Product Type: Dental Laboratory Service

Shipping Options: DHL, UPS, FEDEX

Customer Service: Available 7\*24H

Shade Options: A1, A2, A3, A3.5, and 16 other shades

Flexural Strength: 360 MPa

High esthetic PFZ crowns and bridges for superior appearance

Provides high esthetic dental restorations ensuring natural look and durability

## Technical Parameters:

Service Type	China Dental Laboratory Services
Flexural Strength	360 MPa
Precision Of Margins	Accurate
Resistance To Chipping	High
Strength And Durability	High
Biocompatibility	Safe And Non-toxic Materials
Shipping	DHL/UPS/FEDEX
Feature	Hypoallergenic & Biocompatible
Appearance	Natural-looking
Product Type	Dental Laboratory Service

## Applications:

The ADL Crown And Bridge Dental Lab, originating from China, offers top-tier dental restoration solutions that cater to a wide range of clinical applications. Certified by ISO, CE, and FDA, this product ensures safety, quality, and compliance with international standards, making it a reliable choice for dental professionals worldwide. Its hypoallergenic and biocompatible materials guarantee safe and non-toxic usage, providing peace of mind for both dentists and patients.

One of the primary application occasions for the ADL Crown And Bridge Dental Lab product is in the fabrication of porcelain fused to zirconia crowns and bridges. This technique combines the strength of zirconia with the aesthetic appeal of porcelain, resulting in restorations that are both durable and visually natural-looking. Dental practitioners can confidently utilize this product for cases requiring high strength and superior aesthetics, such as anterior and posterior restorations.

Layered zirconia crowns and bridges are another significant application scenario for this product. The layering process enhances the natural translucency and depth of the restoration, producing a highly aesthetic and natural appearance that closely mimics real teeth. This makes the ADL product ideal for cosmetic dentistry cases where patient satisfaction with the appearance is paramount. Whether it's for single-unit crowns or multi-unit bridges, the layered zirconia technique provides excellent adaptability to various clinical situations.

The ADL Crown And Bridge Dental Lab product is well-suited for dental laboratories and clinics seeking efficient turnaround times, with a delivery time of just 72 hours. The packaging, consisting of secure plastic bags and boxes, ensures that the restorations arrive in pristine condition. Additionally, the service supports flexible payment terms, including Paypal and T/T, and negotiable minimum order quantities, accommodating the needs of both small and large-scale operations.

With a supply ability of 630 technicians, ADL guarantees consistent quality and timely service. This China-based dental laboratory service is dedicated to providing restorations that combine functionality with superior aesthetics, making it an excellent choice for dental professionals aiming to deliver natural-looking and biocompatible dental prostheses. The ADL Crown And Bridge Dental Lab product is truly a comprehensive solution for modern dental restoration needs.

## Support and Services:

Our Crown and Bridge Dental Lab offers comprehensive technical support and services to ensure the highest quality restorations for your

patients. We utilize advanced materials and cutting-edge technology to deliver durable and aesthetically pleasing crowns and bridges tailored to your specifications.

Our team of experienced dental technicians is committed to providing precise craftsmanship and timely delivery. We support a wide range of materials including porcelain-fused-to-metal, all-ceramic, zirconia, and metal alloys to meet diverse clinical needs.

Technical support includes guidance on impression techniques, shade selection, and design considerations to optimize fit and function.

We also offer customized solutions and consultations to address complex cases and special requirements.

To maintain the integrity and performance of your restorations, we provide detailed instructions on handling, cementation, and maintenance. Our customer service team is dedicated to assisting you throughout the process, ensuring satisfaction and successful clinical outcomes.

### Packing and Shipping:

Each Crown and Bridge dental restoration is carefully packaged to ensure maximum protection during transit. The product is placed in a custom-fitted protective case, surrounded by cushioning materials to prevent any movement or damage.

All packages are sealed securely and labeled clearly with handling instructions to maintain product integrity. We use reliable shipping services that offer tracking and timely delivery to ensure your order arrives safely and on schedule.

For international shipments, appropriate customs documentation is included to facilitate smooth clearance. Our packaging and shipping procedures adhere to industry standards to guarantee that your Crown and Bridge dental lab products reach you in optimal condition.



Advanced  
Dental Laboratory



Advanced  
Dental Laboratory

**Advanced Dental Laboratory**



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong  
Province, P.R.China