



China Dental Lab Accurate Precision Of Margins Layered Zirconia Hypoallergenic Biocompatible Dental Material For Durable And Restorations

Our Product Introduction

for more products please visit us on adidentalab.com


Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Resistance To Chipping: High
- Shade: A1/A2/A3/A3.5...16 Shades
- Flexural Strength: 360 MPa
- Strength And Durability: High
- Fit: Precise
- Precision Of Margins: Accurate
- Experience: More Than 13 Years
- Hypoallergenic: Yes



More Images



Product Description

Product Description:

Our F

PFZ is an exceptional solution designed to meet the highest standards of dental restoration. This product is meticulously crafted to provide a natural-looking appearance, ensuring that every crown and bridge blends seamlessly with the patient's existing teeth. Utilizing advanced materials such as layered zirconia and porcelain fused to zirconia (PFZ), the dental lab guarantees restorations that are both high esthetic and functionally superior.

One of the standout features of this product is its extensive shade range. Available in 16 different shades including A1, A2, A3, and A3.5, PFZ offers a wide spectrum of color options to perfectly match the natural tooth color of each patient. This variety allows dental professionals to provide customized restorations that achieve a flawless, natural look, enhancing patient satisfaction and confidence. Durability is a critical factor in dental restorations, and China Dental Lab-Advanced Dental Lab's product excels in this area. Made with long-lasting and wear-resistant materials, these crowns are engineered to withstand the daily rigors of chewing and biting. The use of layered zirconia provides a strong core that supports the restoration, while the PFZ layering offers an esthetic outer surface that resists chipping and wear. This combination ensures that the restorations maintain their integrity and appearance over time, reducing the need for frequent replacements or repairs.

Shipping and delivery are handled with utmost care and efficiency. The product is shipped via reliable carriers such as DHL, UPS, and FEDEX, ensuring timely and secure delivery to dental clinics and dental labs worldwide. The logistics are designed to maintain the product's condition during transit, so dental professionals receive crowns and bridges that are ready for immediate use without any compromise in quality.

Customer service is a cornerstone of the ADL's experience. The lab provides 7*24H support, meaning dental professionals can access assistance at any time of day or night. This round-the-clock customer service ensures that any questions, concerns, or special requests are addressed promptly and professionally. Whether it's technical support, order tracking, or consultation on shade selection, the team is dedicated to providing comprehensive support that enhances the overall service experience.

In summary, China Dental Lab-Advanced Dental Lab offers a perfect blend of natural aesthetics, robust durability, and exceptional service. Its high esthetic appeal, achieved through advanced techniques such as layered zirconia and PFZ, makes it an ideal choice for dental restorations that require both beauty and strength. The wide range of shades ensures an accurate match to natural teeth, while the long-lasting, wear-resistant materials provide peace of mind for both patients and practitioners. Coupled with reliable shipping options and dedicated 24/7 customer service, this product stands out as a top-tier solution in the dental restoration market.

Features:

Product Name: Layered Zirconia/PFZ
Service Type: China Dental Laboratory Full Services
Fit: Precise fit for optimal comfort and durability
Customer Service: Available 7*24H to assist you anytime
Hypoallergenic: Yes, safe for sensitive patients
Shipping: Fast and reliable shipping via DHL, UPS, and FEDEX
Material Options: PFZ and layered zirconia for high esthetic results
High Esthetic: Designed to provide natural-looking crowns and bridges

Technical Parameters:

Product Name	Layered Zirconia/PFZ
Shipping	DHL/UPS/FEDEX
Fit	Precise
Shade	A1/A2/A3/A3.5...16 Shades
Hypoallergenic	Yes
Resistance To Chipping	High
Aesthetic	Highly Aesthetic And Natural Appearance (High esthetic, PFZ)
Precision Of Margins	Accurate
Appearance	Natural-looking
Customer Service	7*24H

Applications:

China dental lab, the ADL Crown And Bridge Dental Lab product, originating from China and certified with ISO, CE, and FDA standards, is designed to meet the diverse needs of dental professionals seeking high-quality dental laboratory services. This product is especially suitable for dental clinics and prosthodontists requiring reliable and precise crown and bridge restorations that offer both durability and a high esthetic appearance. The porcelain fused to zirconia technology used in this product ensures exceptional strength and resistance to chipping, making it an ideal choice for patients who desire long-lasting dental restorations without compromising on natural aesthetics. In various application occasions such as routine dental restorations, cosmetic dental procedures, and rehabilitative dentistry, the Crown And Bridge Dental Lab product by ADL proves to be highly effective. Its high esthetic and highly natural appearance are perfect for anterior and posterior restorations, providing patients with a smile that looks both beautiful and natural. Dental professionals can confidently use this product to fabricate crowns and bridges that blend seamlessly with the surrounding teeth, enhancing overall patient satisfaction.

Scenarios such as repairing damaged teeth, replacing missing teeth, and restoring oral function benefit greatly from the use of ADL's crown and bridge services. The lab's efficient delivery time of 72 hours ensures that dental practices can provide timely treatment, minimizing patient wait times. With a supply ability supported by 630 skilled technicians, ADL can handle both small and large orders, with minimum order quantities negotiable to fit specific requirements. Packaging details include secure plastic bags and boxes, ensuring the safe transport of dental restorations via reliable shipping methods like DHL, UPS, or FEDEX. Payment terms are flexible, accepting Paypal or T/T transactions, making it convenient for dental professionals globally to access these high esthetic restorations. Overall, the ADL Crown And Bridge Dental Lab product combines advanced porcelain fused to zirconia technology with expert craftsmanship to deliver restorations that are not only resistant to chipping but also exhibit a highly aesthetic and natural appearance, suitable for a wide range of dental restoration occasions and clinical scenarios.

Support and Services:

Our Crown and Bridge Dental Lab product comes with comprehensive technical support and services to ensure optimal performance and customer satisfaction. We provide detailed product documentation, including user manuals and maintenance guides, to help you understand and utilize the product effectively.

Our technical support team is available to assist with troubleshooting, installation guidance, and any operational inquiries. We offer remote support services as well as on-site assistance when necessary, ensuring timely resolution of any issues you may encounter.

Additionally, we provide training sessions and workshops for dental professionals and technicians to enhance their skills in using the Crown and Bridge Dental Lab product. Regular software updates and product enhancements are made available to keep your equipment up to date with the latest technology advancements. We are committed to delivering high-quality customer service and technical support to help you achieve the best results in your dental lab work.

Packing and Shipping:

All Crown and Bridge dental lab products are carefully packaged to ensure maximum protection during transit. Each restoration is securely placed in a custom-fitted container with cushioning materials to prevent any damage. Our packaging materials are designed to maintain the integrity and cleanliness of the dental products, ensuring they arrive in perfect condition. We ship all orders using reliable carriers with tracking options available, providing timely delivery and peace of mind to our clients. For international shipments, we comply with all relevant regulations and provide necessary documentation to facilitate smooth customs clearance. Upon receipt, please inspect the package immediately and contact us if there are any issues or damages.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China