



Hypoallergenic PFZ dental solutions tailored to meet diverse clinical needs dental laboratory services ensuring precision

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Customer Service: 7*24H
- Shade: A1/A2/A3/A3.5...16 Shades
- Biocompatibility: Safe And Non-toxic Materials
- Aesthetic: Highly Aesthetic And Natural Appearance
- Productname: Crown And Bridge Dental Lab
- Appearance: Natural-looking
- Fit: Precise
- Producttype: Dental Laboratory Service



More Images



Product Description

Product Description:

Our F

The Crown And Bridge Dental Lab product is designed to meet the highest standards of dental restoration, offering patients a natural-looking and highly aesthetic solution for missing or damaged teeth. This product is meticulously crafted to provide an exceptional blend of durability, precision, and beauty, making it an ideal choice for dental professionals and their patients seeking reliable and visually appealing restorations. Whether you are looking to restore a single crown or multiple bridges, the Crown And Bridge Dental Lab product delivers outstanding results with its advanced technology and expert craftsmanship.

One of the standout features of this product is its natural-looking appearance, which is achieved through the use of porcelain fused to zirconia (PFZ) technology. PFZ combines the strength of zirconia with the exceptional aesthetic qualities of porcelain, resulting in restorations that closely mimic the translucency and color of natural teeth. This ensures that each crown or bridge not only functions effectively but also blends seamlessly with the surrounding dentition, providing patients with a confident and beautiful smile.

The aesthetic appeal of the Crown And Bridge Dental Lab product is highly regarded in the dental community. Thanks to the porcelain fused to zirconia (PFZ) technique, these restorations offer a highly aesthetic and natural appearance that is difficult to distinguish from real teeth. The PFZ framework provides a strong, stable base while the porcelain overlay can be customized with various shades and textures, tailored to match the individual characteristics of each patient's teeth. This level of customization ensures that every restoration is unique and perfectly suited to the patient's oral environment.

In addition to its superior aesthetics, the Crown And Bridge Dental Lab product is known for its precise fit. Precision is critical in dental restorations to ensure comfort, functionality, and long-term success. The lab utilizes state-of-the-art CAD/CAM technology and meticulous quality control processes to guarantee that every crown and bridge fits perfectly within the patient's mouth. This precise fit reduces the risk of complications such as discomfort, gum irritation, or improper bite alignment, allowing for a more natural and comfortable dental experience.

Another significant advantage of the Crown And Bridge Dental Lab product is its reliable and efficient shipping options. To accommodate the needs of dental professionals worldwide, the product is shipped via well-known and trusted courier services including DHL, UPS, and FEDEX. This ensures timely delivery and safe handling of the dental restorations, minimizing downtime and allowing dental practices to provide prompt care to their patients. Fast and secure shipping is essential for maintaining workflow efficiency and patient satisfaction in any dental practice.

Overall, the Crown And Bridge Dental Lab product offers an unparalleled combination of natural-looking aesthetics, precise fit, and durable materials. Its use of porcelain fused to zirconia (PFZ) technology ensures that restorations are both strong and beautiful, meeting the demands of modern dental restoration. Dental professionals can confidently rely on this product to provide their patients with crowns and bridges that not only restore functionality but also enhance the overall appearance of their smiles. Coupled with efficient shipping via DHL, UPS, or FEDEX, the Crown And Bridge Dental Lab product is the perfect choice for dental laboratories and clinics seeking excellence in restorative dentistry.

Features:

Product Name: Crown And Bridge Dental Lab

Biocompatibility: Safe and Non-toxic Materials

Strength and Durability: High

Hypoallergenic: Yes

Aesthetic: Highly Aesthetic and Natural Appearance

Material Options: Porcelain fused to zirconia and layered zirconia

Advanced layered zirconia technology for enhanced strength and esthetics

Technical Parameters:

Resistance To Chipping	High
Flexural Strength	360 MPa
Shade	A1/A2/A3/A3.5...16 Shades
Hypoallergenic	Yes
Chemical Solubility	< 100 µg/mm ²
Customer Service	7*24H
Shipping	DHL/UPS/FEDEX
Experience	More Than 13 Years
Service Type	China Dental Laboratory Services
Product Type	Dental Laboratory Service

Applications:

The ADL Crown And Bridge Dental Lab product is an ideal solution for a variety of dental restoration applications, offering exceptional quality and performance. Manufactured in China, this product is certified by ISO, CE, and FDA, ensuring compliance with international safety and quality standards. The lab specializes in creating high-precision dental restorations such as layered zirconia crowns and bridges, as well as porcelain fused to zirconia (PFZ) restorations, which are renowned for their high esthetic appeal and durability.

This product is widely used in dental clinics and laboratories for patients requiring natural-looking prosthetics that blend seamlessly with their existing teeth. The crowns and bridges fabricated by ADL feature accurate margin precision, which is crucial for achieving a perfect fit and preventing bacterial infiltration. With a flexural strength of 360 MPa, these restorations offer excellent resistance to chipping and mechanical stress, making them suitable for both anterior and posterior applications.

ADL's Crown And Bridge Dental Lab product is highly valued in scenarios where biocompatibility is paramount. Made from safe and non-toxic materials, these restorations ensure patient safety and comfort over long-term use. The high esthetic qualities of layered zirconia and PFZ restorations provide a natural-looking appearance, making them an excellent choice for patients seeking not only functional but also visually appealing dental solutions.

The product is conveniently packaged in plastic bags and boxes to maintain hygiene and facilitate easy handling. ADL offers flexible purchasing options with negotiable minimum order quantities and prices, catering to both small-scale and large-scale dental labs. The lab boasts a supply ability capable of supporting 630 technicians, ensuring timely delivery within 72 hours. Payment can be securely made via Paypal or T/T, providing convenience and reliability for customers worldwide.

Overall, the ADL Crown And Bridge Dental Lab product is perfect for dental professionals looking for advanced, high-quality restorations that combine strength, precision, and esthetics. Whether for standard crowns, multi-unit bridges, or specialized PFZ restorations, this product meets the demands of modern dental practice with outstanding results.

Support and Services:

Our Crown and Bridge Dental Lab product is designed to provide dental professionals with high-quality, precise, and durable restorations. For technical support, we offer comprehensive assistance to ensure optimal use and performance of our product. Our services include troubleshooting guidance, detailed product usage instructions, and maintenance recommendations to maximize the lifespan and effectiveness of your dental restorations.

We also provide training resources and expert consultations to help dental technicians and clinicians achieve the best results with our

Crown and Bridge solutions. Whether you need help with material selection, digital workflow integration, or finishing techniques, our team is dedicated to supporting your success.

Additionally, our technical support extends to software updates and compatibility checks to keep your digital dental lab operations running smoothly. We encourage users to utilize our online knowledge base and FAQs for quick answers and practical tips.

Our commitment is to deliver exceptional customer service and ensure that your experience with our Crown and Bridge Dental Lab product is seamless and productive. For further assistance, please refer to our support resources and service options.

Packing and Shipping:

All Crown and Bridge Dental Lab products are carefully packaged to ensure maximum protection during transit. Each restoration is securely placed in a custom-molded protective case, surrounded by cushioning materials to prevent movement and damage.

Our packaging materials are designed to maintain the integrity and cleanliness of the dental restorations, keeping them free from contamination.

Once packaged, the products are sealed in a sturdy, tamper-evident outer box labeled with handling instructions and product details.

We offer reliable shipping options with tracking to ensure timely and safe delivery to your dental practice or laboratory. Expedited shipping is available upon request to meet urgent needs.

Please inspect the package immediately upon arrival and notify us of any damage or discrepancies so we can address the issue promptly.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China