



Innovative Dental Lab Outsourcing Services Combining Advanced Techniques And Control For Dental Healthcare Providers

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information



Product Specification

- Highlight: dental lab outsourcing services with warranty, China dental lab advanced techniques, dental healthcare provider outsourcing solutions



Our Pr

Product Description

Dental Lab Outsourcing Services for All-on-X Implant Retained Bar | Screw Retained Zirconia Bridge

Product Description

Our **dental lab outsourcing services** specialize in high-performance All-on-X implant retained bar solutions combined with screw-retained zirconia bridges. Engineered for superior retention and long-term stability, this system is ideal for full-arch restorations requiring strength, precision, and durability.

Through advanced CAD/CAM workflows, our **dental lab outsourcing services** ensure precise fabrication of both the titanium bar and zirconia bridge, delivering passive fit and minimizing clinical adjustments. Each case is designed to meet the demanding requirements of modern implant dentistry.

By choosing professional **dental lab outsourcing services**, dental clinics and laboratories can enhance efficiency, reduce in-house workload, and maintain consistent quality across complex implant restorations.

Key Features & Advantages

Feature	Description
Restoration Type	All-on-X implant retained bar with screw-retained zirconia bridge
Retention & Stability	Engineered for superior retention and long-term structural stability
Material Combination	Titanium bar + full zirconia bridge for strength and aesthetics
Manufacturing Technology	Precision CAD/CAM production through reliable dental lab outsourcing services
Fit Accuracy	Passive fit design reduces stress on implants and minimizes adjustments
Aesthetic Performance	Natural-looking zirconia with customized shading options
Customization	Tailored solutions based on individual case requirements
Quality Control	Multi-stage inspection ensures consistent and reliable results
Production Efficiency	Streamlined workflow supported by professional dental lab outsourcing services
Cost Advantage	Competitive pricing with scalable outsourcing solutions

Engineered for Implant Success

This system combines a rigid implant-retained bar with a screw-retained zirconia bridge, offering enhanced mechanical stability and optimal load distribution. Designed by experienced teams using advanced **dental lab outsourcing services**, it ensures reliable outcomes even in complex full-arch cases.

Ideal for Modern Dental Practices

Our **dental lab outsourcing services** are tailored for clinics and laboratories seeking high-quality implant restorations with predictable performance. The combination of strength, precision, and efficiency helps reduce chair time and improve patient satisfaction.

Applications

- Full-arch All-on-4 / All-on-6 restorations
- Implant-retained bar systems
- Screw-retained zirconia bridges
- Complex implant-supported prosthetics

All restorations are produced under strict standards using advanced **dental lab outsourcing services**, ensuring long-term clinical reliability.



Advanced Dental Laboratory



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China