



Professional Dental Lab Outsourcing Services Delivering Accurate Dental Prosthetics and Custom Solutions for Dental Industry Needs

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information



Product Specification

- Highlight: professional dental lab outsourcing services, accurate dental prosthetics solutions, custom dental industry prosthetics



Product Description

Our Pr

China Outsourcing Dental Lab Titanium Full Arch Hybrid Bridge | Immediate Loading A1 Shade Screw Retained CAD/CAM

Product Description

As a professional **china outsourcing dental lab**, we specialize in manufacturing titanium full arch hybrid bridges designed for immediate loading protocols. Combining advanced CAD/CAM technology with precise engineering, our restorations ensure superior strength, accurate fit, and long-term clinical reliability.

Our **china outsourcing dental lab** produces screw-retained hybrid bridges featuring a durable titanium framework with aesthetic zirconia or composite layering. Designed in A1 shade with refined Malong-style morphology, these restorations deliver a natural appearance while maintaining structural integrity.

By working with an experienced **china outsourcing dental lab**, dental clinics and laboratories can reduce production costs, improve turnaround efficiency, and maintain consistent quality across full-arch restorations.

Key Features & Advantages

Feature	Description
Restoration Type	Titanium full arch hybrid bridge (screw-retained)
Loading Protocol	Immediate loading supported for faster treatment workflows
Material Structure	Titanium framework with zirconia/composite layering
Shade & Design	A1 shade with natural Malong-style anatomical design
Manufacturing Technology	Precision CAD/CAM production from a reliable china outsourcing dental lab
Fit Accuracy	Passive fit reduces implant stress and chairside adjustments
Customization	Tailored designs based on individual clinical requirements
Quality Control	Multi-stage inspection ensures consistent performance
Production Efficiency	Streamlined workflow supported by a professional china outsourcing dental lab
Cost Advantage	Competitive factory pricing with scalable outsourcing solutions

Designed for Immediate Loading

This titanium full arch hybrid bridge is engineered for immediate loading cases, providing strong primary stability and predictable outcomes. Developed by a trusted **china outsourcing dental lab**, it helps clinicians deliver faster and more efficient treatments.

Aesthetic & Functional Performance

With A1 shade selection and detailed anatomical contouring, the restoration achieves a highly natural look. Manufactured using advanced CAD/CAM systems in a certified **china outsourcing dental lab**, each case meets high standards for both aesthetics and function.

Applications

- Full-arch implant restorations
- Immediate loading cases
- Screw-retained hybrid bridges
- All-on-4 / All-on-6 prosthetic solutions

All restorations are produced under strict standards by an experienced **china outsourcing dental lab**, ensuring reliable long-term clinical performance.



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China