

Medical-grade Nylon Occlusal Splint Night Guard with Flexible CAD/CAM Design for Bruxism Protection

Our Product Introduction

for more products please visit us on addentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Model Number: ND-01
- Minimum Order Quantity: 1 Piece
- Packaging Details: Individual sealed box with protective case, export-grade carton
- Delivery Time: 7 days
- Payment Terms: T/T, PayPal



Product Specification

- Material: Medical-grade Nylon (thermoplastic)
- Splint Type: Dorsal Occlusal Splint
- Hardness: Flexible / Semi-rigid
- Thickness: 1.5-2.5 Mm
- Fabrication: CAD/CAM Digital Design And Milling
- Color: White
- Suitable For: Bruxism, TMJ Disorders, Teeth Grinding Protection
- Highlight: **Flexible Nylon Dorsal Splint, Bruxism Night Guard Protection, Custom Occlusal Splint Lab**



More Images



Product Description

Nylon Dorsal Flexible Occlusal Splint Bruxism Night Guard for Teeth Grinding Protection

Our P

Product Description:

The Nylon Dorsal is a flexible occlusal splint fabricated from medical-grade nylon, designed to protect the teeth from grinding and clenching during sleep. Its dorsal design covers the occlusal surface of the dentition, providing a comfortable barrier that absorbs excessive bite force and relieves pressure on the temporomandibular joint.

Unlike bulky rigid acrylic guards, the Nylon Dorsal is lightweight and flexible, allowing it to fit snugly and comfortably in the mouth. Each splint is custom-made from a dental impression or intraoral scan, so it stays securely in place throughout the night without slipping or causing discomfort.

Thanks to its durable nylon material, the splint offers excellent long-term performance. It resists wear and maintains its shape even with daily use, making it a reliable solution for patients who suffer from bruxism and jaw muscle tension.

Overall, the Nylon Dorsal is an excellent choice for anyone looking for a comfortable, precise-fitting occlusal splint. Its flexible nylon construction, custom fit, and clean white appearance make it a discreet and effective way to protect natural teeth and restorations from the harmful effects of teeth grinding.

Features:

Product Name: Nylon Dorsal Occlusal Splint

Material: Medical-grade nylon

Color: White

Splint Type: Dorsal occlusal splint

Fit: Custom-made from impression or intraoral scan

This splint is designed to protect the teeth and jaw from the effects of bruxism. Its dorsal occlusal coverage helps to distribute bite force evenly, while the flexible nylon material ensures a comfortable and secure fit for all-day or overnight wear.

Technical Parameters:

Material	Medical-grade nylon (thermoplastic)
Splint Type	Dorsal occlusal splint
Hardness	Flexible / semi-rigid
Thickness	1.5-2.5 mm
Fabrication	CAD/CAM digital design and milling
Color	White

Applications:

The Nylon Dorsal is ideal for patients who grind or clench their teeth while sleeping, as well as those with temporomandibular joint disorders and jaw muscle tension. It is also suitable for protecting dental restorations and natural dentition from excessive wear.

The Nylon Dorsal is manufactured in Shenzhen, China, under strict quality control measures. It is certified by ISO, CE, and FDA, which guarantees its safety and quality. The product is individually packed in a protective case to ensure it remains clean and secure during transportation and storage.

ADL has a strong supply capacity to meet both single and bulk orders, and the minimum order quantity is negotiable. Customers can pay using T/T or PayPal, and the delivery time is 7 days, making it convenient for clinics and laboratories that need a reliable, fast turnaround. In conclusion, the ADL Nylon Dorsal is a high-quality occlusal splint that offers excellent protection for patients who grind their teeth. Its flexible nylon material, custom fit, and clean white appearance make it comfortable to wear, while its certifications and fast delivery make it a reliable choice for dental professionals.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China